


Infusomat[®]plus Line Portfolio

Product Information


B

BRAUN

SHARING EXPERTISE

Infusomat [®] plus Line portfolio	Description	Tube material	Length (cm)	Priming vol. (ml)	Y-port	Code no. (REF)
Enteral Nutrition – ENFit™ Lines						
	▪ Inline drip chamber ▪ Different tube materials ▪ Different lengths ▪ ENFit™ patient connector	PUR**	330/220	28.0	needle-free	8700370
	▪ ENPlus Spike for enteral nutrition diets with ENPlus (cross) port	PVC*	330/220	23.1	needle-free	8700400
	▪ Stopcock with male ENFit™ connector for enteral syringe with female ENFit™ connector	PVC*	260/150	18.5	needle-free	8700380
Regional Anesthesia – NRFit™ Line						
	▪ 15 µm fluid filter ▪ Drip chamber incl. AirStop function ▪ NRFit™ patient connector ▪ translucent yellow tube	PUR**	300/200	21.8	–	8700410

* PVC = DEHP-free; ** PUR = PVC-/DEHP-free



Infusomat[®]plus Line Portfolio

Product Information

Infusomat[®] plus Line Portfolio

Product Information

Time to Care

The color codes and shapes of Infusomat[®] compact^{plus} pump make it easy to insert the Infusomat[®] plus lines. The dedicated tubing material allows continued accuracy over long periods (96 h). The Anti Free-Flow Clamp automatically closes the line to prevent any risk of free-flow when the door is opened. In addition, the graphical user interface allows intuitive operation and guided loading reduces training needs.

Infusomat [®] plus Line portfolio	Description	Tube material	Length (cm)	Priming vol. (ml)	Y-port	Code no. (REF)
SafeSet Lines						
	15 µm fluid filter	PVC*	240/150	17.3	–	8700200
	Drip chamber incl. AirStop function	PVC*	300/200	21.2	–	8700210
	Different tube materials	PUR**	240/150	17.8	–	8700220
	Different lengths	PVC*	300/200	22.2	needle-free	8700240
	Luer-Lock Patient connector	PUR**	300/200	22.9	needle-free	8700250
Standard Lines						
	15 µm fluid filter	PVC*	240/150	17.3	–	8700310
	Spike drip chamber	PUR**	240/150	17.8	–	8700320
	Different tube materials	PVC*	300/200	22.2	needle-free	8700330
	Luer-Lock patient connector	PVC*	300/200	21.7	needle-based	8700340
Transfusion Lines						
	200 µm blood filter	PVC*	240/150	17.3	–	8700350
	Spike drip chamber					
	Different lengths					
	Spinlock patient connector	PVC*	300/200	21.9	needle-free	8700360
UV-light protect Lines						
	15 µm fluid filter	PUR**	240/150	17.8	–	8700260
	Drip chamber incl. AirStop function					
	Different lengths	PUR**	300/200	22.9	needle-free	8700270
	Luer-Lock patient connector					
Oncology Line						
	15 µm fluid filter					
	0.2 µm hydrophilic fluid filter (Sterifix [®])					
	Drip chamber incl. AirStop function	PUR**	240/150	20.6	–	8700300
	Eurofix Spinlock patient connector					

* PVC = DEHP-free; ** PUR = PVC-/DEHP-free

Infusomat [®] plus Line portfolio	Description	Tube material	Length (cm)	Priming vol. (ml)	Y-port	Code no. (REF)
Cyto-Set [®] Lines						
	Standard Lines:					
	15 µm fluid filter	PUR**	210/155	23.9	3 needle-free valves	8700480
	Grip plate with integrated back-check valves					
	Cytostatic resistant material	PUR**	210/155	23.9	5 needle-free valves	8700490
	Eurofix Spinlock patient connector					
	0.2 µm hydrophilic fluid filter (Sterifix [®])	PUR**	225/170	27.4	5 needle-free valves	8700470
	AirStop Lines:					
	Drip chamber incl. AirStop function	PUR**	210/155	23.9	3 needle-free valves	8700420
	Grip plate with integrated back-check valves					
	Cytostatic resistant material	PUR**	210/155	23.9	5 needle-free valves	8700430
	Eurofix Spinlock patient connector					
	0.2 µm hydrophilic fluid filter (Sterifix [®])	PUR**	225/170	27.4	5 needle-free valves	8700460
	UV-light protect Lines:					
	Light protection up to 520 nm	PUR**	210/155	23.9	3 needle-free valves	8700440
	Drip chamber incl. AirStop function					
	Grip plate with integrated back-check valves					
	Cytostatic resistant material	PUR**	210/155	23.9	5 needle-free valves	8700450
	Eurofix Spinlock patient connector					
Flushing Set						
	15 µm fluid filter					
	Inline drip chamber incl. AirStop function	PUR**	300/200	22.4	needle-free	8700280
	Spinlock patient connector					

* PVC = DEHP-free; ** PUR = PVC-/DEHP-free