


Original KLEINSASSER Operating Laryngoscopes

Full size illustration:
inner diameter in mm proximal and distal



8590 A - L

8590 GF - HF

8590 A KLEINSASSER Operating Laryngoscope, for adults, extra large, length 17 cm

8590 B Same, large

8590 C Same, medium (most commonly used model)

8590 DN Same, small, length 18 cm (for particularly difficult anatomical circumstances)

8590 J Same, medium, length 18 cm (for anterior commissure)

8590 JL Same, medium, length 22 cm (for especially long neck)

8590 K Same, for children, small, length 15 cm

8590 L Same, for infants, small, length 13 cm

497 AC Clip, for proximal illumination


8590 GF Fiber Optic Light Carrier, for distal illumination, working length 14 cm, for use with Laryngoscopes 8590 A - JL

8590 HF Fiber Optic Light Carrier, for distal illumination, working length 7.5 cm, for use with Laryngoscopes 8590 K - L

8590 KC Injection Cannula, for positive pressure assisted ventilation, working length 8.5 cm, for use with Laryngoscopes 8590 A - L

8590 M Suction Tube to remove vapor, for LASER treatment, diameter 3 mm, working length 12.5 cm, for use with Laryngoscopes 8590 A - JL

8590 N Suction Tube, to remove vapor, for LASER treatment, diameter 3 mm, working length 9 cm, for use with Laryngoscopes 8590 K - L



Lot 30