

Lot 2.5 si 2.6

Endo GIA™ reloads with Tri-Staple™ technology

Endo GIA™ curved tip reload with Tri-Staple™ technology



Order Code	Color Code	Description	Box Qty
Lot 2.5 EGIA30CTAV	■	Endo GIA™ curved tip 30mm reload articulating vascular reload*	6
EGIA30CTAVM	■	Endo GIA™ curved tip 30 mm articulating vascular/medium reloads With Tri-Staple™ technology	6
Lot 2.6 EGIA45CTAV	■	Endo GIA™ curved tip 45 mm articulating vascular reload*	6
EGIA45CTAVM	■	Endo GIA™ curved tip 45 mm articulating vascular/medium reloads With Tri-Staple™ technology	6
EGIA60CTAVM	■	Endo GIA™ curved tip 60 mm articulating vascular/medium reloads With Tri-Staple™ technology	6
EGIA60CTAMT	■	Endo GIA™ curved tip 60 mm articulating medium/thick reloads With Tri-Staple™ technology	6

* Code EGIA30CTAV and EGIA45CTAV vascular reloads do not come with Tri-Staple™ technology, have a slimmer fixed anvil