


Original KLEINSASSER Operating Laryngoscopes

Full size illustration:
inner diameter in mm proximal and distal




8590 A


8590 B




8590 C



8590 J - JL



8590 DN




8590 K



8590 L



497 AC



8590 GF - HF



8590 KC



8590 M - N



8590 A - L

8590 GF - HF

8590 A **KLEINSASSER Operating Laryngoscope**, for adults, extra large, length 17 cm

8590 B **Same**, large

8590 C **Same**, medium (most commonly used model)

Lot 43

8590 DN **Same**, small, length 18 cm (for particularly difficult anatomical circumstances)

8590 J **Same**, medium, length 18 cm (for anterior commissure)

8590 JL **Same**, medium, length 22 cm (for especially long neck)

8590 K **Same**, for children, small, length 15 cm

8590 L **Same**, for infants, small, length 13 cm

497 AC **Clip**, for proximal illumination

8590 GF **Fiber Optic Light Carrier**, for distal illumination, working length 14 cm, for use with Laryngoscopes 8590 A - JL

8590 HF **Fiber Optic Light Carrier**, for distal illumination, working length 7.5 cm, for use with Laryngoscopes 8590 K - L

8590 KC **Injection Cannula**, for positive pressure assisted ventilation, working length 8.5 cm, for use with Laryngoscopes 8590 A - L

8590 M **Suction Tube to remove vapor**, for LASER treatment, diameter 3 mm, working length 12.5 cm, for use with Laryngoscopes 8590 A - JL

8590 N **Suction Tube, to remove vapor**, for LASER treatment, diameter 3 mm, working length 9 cm, for use with Laryngoscopes 8590 K - L