

MX-S+ Digital Vortex Mixer

Features

- LCD screen displays real-time speed and time for easy adjustment and monitoring.
- Adjustable speed ranges from gentle to vigorous mixing to suit various needs.
- Stable with pressure-mounted accessories to prevent slippage and splashing.
- Brushless DC motor for high stability and maintenance-free operation.
- Sturdy base ensures stability during operation for safety.
- Includes various accessories for different sizes and types of tubes and plates.
- Low noise for a quiet working environment.





Specifications	MX-S+
Speed	100-3000 rpm
Display	LCD
Amplitude	4 mm
Operating Modes	Touch, Continuous
Time Range	1 s - 99 min 59 s
Motor Type	Brushless DC Motor
Dimensions (L \times W \times H)	160x146x115mm
Weight	3.4kg
Power Supply	100-240 V, 50/60 Hz
Input Voltage	24 V
Power Consumption	25 W
Ambient Temperature & Humidity	5-40°C, 80% RH
Enclosure Protection	IP21

Accessories of MX-S+&MX-Pro

VT2.3

Cat. No.18901365

Universal top plate, compatible with VT1.3.1-VT1.3.12 Tube Adapter



VT2.5

Cat. No.18901367

Erlenmeyer flask attachment (100mL-250mL), Operating speed: 0~1500rpm



VT2.4

Cat. No.18901366

Ampule/Tube Attachment, 4pc ampule/tube (0.2mL-15mL)



PS2.1

Cat. No.18901360

Microplate clamp, Suitable for workboards (full skirt PCR plates, cell culture plates, enzyme-linked immunosorbent assay plates, deep well plates, etc. Half skirt/no skirt PCR plates need to be used with rectangular foam, which is already included)



Accessories of MX-S+&MX-Pro (New)

VT2.1

Cat. No.18901379

Top cup, Suitable for MX-Pro



VT2.6.4

Cat. No.18901375

Tube Adapter,8 holes ×15mL Test Tube



VT2.1

Cat. No.18901380

Top cup, Suitable for MX-S+



VT2.6.5

Cat. No.18901374

Tube Adapter, 4 holes × 50mL Test Tube



VT2.2

Cat. No.18901373

Flat Rubber Sleeve (Standard)



VT 2.7.1

Cat. No.18901362

Horizontal Plastic Clamp, 24×1.5mL Test Tube



VT2.6

Cat. No.18901368

Snap-in Universal top plate, compatible with VT1.3.1-VT1.3.12 Tube Adapter



VT 2.7.2

Cat. No.18901363*

Horizontal Plastic Clamp, 12×15mL Test Tube



VT2.6.1

Cat. No.18901378

Tube Adapter, $25 \text{ holes} \times 1 \text{mL}$ Test Tube



VT 2.7.3

Cat. No.18901364*

Horizontal Plastic Clamp, 6×50mL Test Tube



VT2.6.2

Cat. No.18901377

Tube Adapter, 21 holes $\times 2 \mathrm{mL}$ Test Tube



PS2.2

Cat. No.18901361

Snap-in Microplate clamp, Suitable for workboards (full skirt PCR plates, cell culture plates, enzyme-linked immunosorbent assay plates, deep well plates, etc. Half skirt/no skirt PCR plates need to be used with rectangular foam, which is already included)



VT2.6.3

Cat. No.18901376

Tube Adapter, 16 holes $\times 5 \mathrm{mL}$ Test Tube



^{*}Material color may vary slightly during the initial launch phase. Please refer to the actual product.