

HANGZHOU SHANYOU MEDICAL EQUIPMENT CO., LTD
No. 138, Louta Development Zone, Guancun Village, Louta
Town,
Xiaoshan District, Hangzhou, Zhejiang 311266, China
Tel: 86-571-82304593 Fax: 86-571-82305195

Standard Endotracheal Tube


The endotracheal tube is inserted into the trachea through glottis to establish and maintain an artificial airway when take the anesthesia and emergency treatment.



Features:

1. Made of clear, soft, medical-grade thermosensitive PVC, high quality,
2. Marking line in cm
3. High-volume, low-pressure cuff, with Murphy eye
4. Radiopaque line throughout the tube for X-ray Visualization, Transparent
5. Sterile, EO
6. Individually packaged
7. CE, ISO Certified
8. Maximum pressure 40mmHg
9. Well placed connector
10. Good contact with trachea

Standard Endotracheal Tube	
Cuffed	
Size	Product No.
2.5	ET1P25C
3	ET1P30C
3.5	ET1P35C
4	ET1P40C
4.5	ET1P45C
5	ET1P50C
5.5	ET1P55C
6	ET1P60C
6.5	ET1P65C
7	ET1P70C
7.5	ET1P75C
8	ET1P80C
8.5	ET1P85C
9	ET1P90C
9.5	ET1P95C
10	ET1P10C



李石明