

Test report No. 24021sd

EVALUATION OF FUNGICIDAL OR YEASTICIDAL ACTIVITY IN THE MEDICAL AREA (EN 13624)

Name of the product*:

BACTICID AF

Batch number*:

197271120

Order number:

20037

Manufacturer*:

Chemi-Pharm Ltd.

Client, representative*:

Chemi-Pharm Ltd., Tänassilma tee 11, Tänassilma küla, Saku vald,

76406, ESTONIA, Maris Millner, +372 5177090

Date of delivery:

10.12.2020

Test material conditions:

No specific features, sample in the manufacturers tare

Storage conditions:

At room temperature, darkness

Active substance – conc.*:

Ethyl alcohol 57.0% wt, isopropyl alcohol 6.0% wt

Appearance of the product:

Transparent liquid

Test concentration:

80.0%, 50.0%, 25.0%

Contact time:

2 min

Interfering substance:

3g/I bovine albumin + 3 ml/I sheep blood erythrocytes (dirty

conditions)

Neutralizer:

Rinsing liquid:

Tryptone 1 g/l + NaCl 9 g/l

Test organisms:

Aspergillus brasiliensis ATCC 16404

Testing method:

EVS-EN 13624:2013

Quantitative suspension test for the evaluation of fungicidal or

yeasticidal activity in the medical area.

Testing date:

17.02.2021 - 19.02.2021

Results:

Look appendix 1

Interpretation and conclusion: Look appendix 2

Allar Laaneleht Chief specialist

Date of issue: 26.04.2021

* - Data provided by the customer



E-mail: info@ntl.ee



Test report No. 24021sd

Appendix 1

TEST RESULTS (suspension test)

EVS-EN 13624:2013; Phase 2, step 1

Membrane filtration method

Rinsing liquid: Tryptone 1 g/l + NaCl 9 g/l;

Test organism: Aspergillus brasiliensis ATCC 16404

Test temperature: $+20^{\circ}$ C; Incubation temperature: $+30 \pm 1^{\circ}$ C

Interfering substance: 3g/l bovine albumin + 3 ml/l sheep blood erythrocytes

Nordic Tersus Laboratory LLC.; Date of test: 17.02.2021

Responsible person: Melissa Ingela Bramanis, Allar Laaneleht

Validation and controls

Dirty conditions

Validation suspension N _{vo}			Experimental conditions (A)			Filtration control (B)			Method validation (C)		
V _{C1}	V _{C2}	Σ̈́	V _{C1}	V _{C2}	Ϋ́	V _{C1}	V _{C2}	Ā	V _{C1}	V _{C2}	x
43	40	41.5	36	30	33	39	42	40.5	45	49	47
30 ≤ x̄ N _{vo} ≤160?yes X; no □			$\bar{x} \mathbf{A} \text{ is } \geq 0.5 \bar{x} \mathbf{N}_{vo} \text{?yes X;no } \square$			\bar{x} B is $\geq 0.5 \bar{x}$ N_{vo} ? yesX; no \Box			x C is ≥ 0.5 x N _{vo} ? yes X;no □		

Test suspension and test

Testsuspension:	N	V _{C1}	V_{C2}	$\bar{x}_{wm} = 2.64 \times 10^7$; $\log N = 7.42$
	10 ⁻⁵	247	275	$N_0 = N/10$; $\log N_0 = 6.42$
N and N_0	10-6	28	28	6.17≤ log N ₀≤6.70; yesX; no □

Experimental results

Concentration of the product %	Dilution step	V _{C1}	V _{C2}	Na (=x̄*10)	lg <i>Na</i>	<i>lg</i> R	Contact time	Conditions
80.0%	-	19	21	200	2.30	4.12	2 min	Dirty
50.0%	-	>55	>55	>550	>2.74	<3.54	2 min	Dirty
25.0%	-	>55	>55	>550	>2.74	<3.54	2 min	Dirty

Explanations:

 V_C = count per ml (one plate or more)

 \bar{x} = average of V_{C1} and V_{C2} (1. + 2. Duplicate)

N = cfu/ml microbes in test suspension

 N_0 = cfu/ml at the start of the contact time (t=0)

 N_{vo} = cfu/ml in the validation suspension (t=0)

Na = surviving microbes after the test

R = reduction factor (R= N_0/Na ; LogR=Log N_0 - Log Na)



Appendix 2

Interpretation:

The ready to use surface disinfection product **BACTICID AF** (batch no. 197271120) was tested according to the test method EVS-EN 13624:2013. The test was performed at 20 °C \pm 1 °C, under dirty conditions during contact times of 2 min. The membrane filtration method was used for testing the product's effectiveness against the reference strains *Aspergillus brasiliensis* ATCC 16404. Under dirty conditions, the 80.0% solution of the tested sample of the product was effective against *Aspergillus brasiliensis* within 2 minutes.

Conclusion:

The surviving count of the reference strain showed at least 4 lg reduction meaning that according to EVS-EN 13624:2013 the sample of the ready to use surface disinfection product BACTICID AF is effective against *Aspergillus brasiliensis* within 2 minutes.

This is the end of the test report

Allar Laaneleht Chief specialist

Date of issue: 26.04.2021