



Emboshield NAV⁶


Embolic Protection System

Navigation

Access

Visibility

Navigate
with a New Level of Control



Navigate with a New Level of Control

Emboshield NAV⁶ Embolic Protection System

Navigation

- Distal Access:** A softer BareWire with extra flexibility and torqueability for successful access through tortuous anatomy.
- Workhorse:** Provides optimal flexibility and support for superb access through most anatomies.
- Support:** Designed to provide additional support for delivery of filter in challenging arches.


Access

Three filter delivery options designed to provide a customized approach to access diverse vascular anatomies.


- | | | |
|--|--|---|
| BareWire: Wire access followed by easy filter delivery through tight lesions and challenging anatomies. | Fixed: Preloaded wire filter delivery system for one-step advancement through lesion. | Telescope: Alternate advancement of BareWire and delivery system for outstanding access and support. |
|--|--|---|

Visibility

Radiopaque circumferential frame conveys vessel apposition and deployment status.



Short Filter Length
19.0 and 22.5 mm



NAV⁶

Advancing Carotid Technology

Emboshield NAV⁶

Embolic Protection System

Navigation / Access / Visibility

Ordering Information

Emboshield NAV⁶ Embolic Protection System*

Unconstrained Filter Diameter	Vessel Diameter	Stock Number	Wire Length	Minimum Sheath/Guide Size	Minimum Sheath/Guide I.D. Required
5.0 mm	2.5–4.8 mm	22442-19 (small)	190 cm	5F / 6F	0.070"
7.2 mm	4.0–7.0 mm	22443-19 (large)	190 cm	5F / 6F	0.070"


*Emboshield NAV⁶ Embolic Protection System consists of 2 primary packs:

Primary Pack	Content
Filtration System	(1) RX Delivery Catheter, (1) Filtration Element, (1) BareWire Filter Delivery Wire (Workhorse, 190 cm), (1) Loading Funnel, (1) Introducer, (1) Torque Device, (1) Syringe with Flushing Tip
RX Retrieval Catheter	(1) RX Retrieval Catheter, (1) Torque Device, (1) Syringe Flushing Tip

BareWire Filter Delivery Wires (0.014")†

Wire Type	Wire Length	Stock Number
Distal Access	190 cm	22444-19
Workhorse	190 cm	22445-19
	315 cm	22445-31
Support	190 cm	22446-19

†BareWires sold in packs of five.



Please contact your local representative for more information.

Abbott Vascular International BVBA

Park Lane, Culliganlaan 2B, 1831 Diegem, Belgium, Tel: 32.2.714.14.11

Tests performed by and data on file at Abbott Vascular. All drawings are artist's representations only and should not be considered as an engineering drawing or photograph. Photo(s) on file at Abbott Vascular. **All products herein** are trademarks of the Abbott Group of Companies.

Caution: This product is intended for use by or under the direction of a physician. Prior to use, it is important to read the package insert thoroughly for Instructions for Use, Warnings and Potential Complications associated with use of this device. Information contained herein for Europe, Middle East and Africa ONLY. Please check the regulatory status of the device before distribution in areas where CE marking is not the regulation in force. For more information, visit our website at abbottvascular.com

