



**One Piece Ostomy Bag**

Stoma care



MODEL	ITEM CODE	DESCRIPTION
Closed Type	1101	Skin Film, Non woven, Without filter, 32mm-60mm
	1102	Skin Film, Non woven, With filter, 32mm-60mm
Drainable Type	1601	Integrated Closure, Transparent Film, Non woven, Without filter, 15mm-60mm
	1611	Integrated Closure, Skin Film, Non woven, Without filter, 15mm-60mm
	1802	Velcro Closure, Transparent Film, Non woven, With filter, 15mm-60mm
	1811	Velcro Closure, Skin Film, Non woven, Without filter, 15mm-60mm
	1512	Clip Type, Transparent Film, Non woven, Without filter, 15mm-60mm
	1381	Pediatric Type(diamond), Velcro Closure, Transparent Film, Non woven, Without filter, 10mm-35mm
Urostomy Type	1010	Transparent Film, Non woven, 10mm-50mm
	1011	Transparent Film, Non woven, 15mm-50mm

**ADVANTAGES**


1. hydrocolloid Flange. It's skin-friendly, hypoallergenic, anti-leaking and prevents edge warping and protects against soaking. The colostomy bag can be stuck on the skin for 3 to 7 days.
2. Non-woven outer lining is provided for the colostomy bag, no super-heating and under cooling feelings, reduce noise when chaling.
3. Filter. It eliminates odor and gas build-up more effectively.

**INTENDED USE**

To collect feces for patients with ileum, rectum and colon operations.

**PACKAGE**

300pcs in a carton (10pcs/box, 30boxes/carton)



Integrated type  
1612

