

With this statement,

as the manufacturer of the patient monitors

we, **Edan Instruments Inc.** do confirm that offered devices comply the following technical specifications:

for the offered monitor IX15

- Capnography (CO₂) Mode, Apnea Recording Time, 10 - 40 s
Criteria for the onset of apnea registration:

In anesthesiology and intensive care, apnea is usually registered if the absence of breathing lasts **≥ 10 seconds**. Complete absence of a capnographic curve (zero EtCO₂) for **≥ 10 seconds**. These are the standards for recording of Apnea in capnography. Critical duration: from 1 minute and more → risk of severe hypoxia, brain damage and cardiac arrest.

- It has Alarm in case of No external power source

Date: 11th, June, 2025



Name: Alvin Liang
Title: Regional Manager