

Quality time for better care

"Quality time for better care" is the brand promise of Terumo Medical Care Solutions.



HOSPITAL CARE

Administration Sets



Supporting your
infusion therapy

A clear view on infusion with a comprehensive portfolio of solutions



Easy handling

- Easy setting following the instruction screen of the pump
- Sharp 2-way spike to penetrate bags and bottles easily
- Adjustable roller clamp with notch to mount tubing during priming

Increased patient safety

- AFF clip for prevention of accidental bolus infusion*
- Does not contain DEHP
- Y-site with NLCD* for prevention of needle stick injuries*

Highly accurate infusion

- Reliable high accuracy combination of Terufusion™ pump & sets[#]
- Durable tubes for continuous accurate infusion[§]

* Needleless connecting device

* Available on specific sets

[#] Flow rate accuracy: within $\pm 5\%$

[§] Always consult the product label and IFU for a complete overview of warnings, cautions and/or precautions prior to actual use.



General

Product code	Ti*PU300WY	Ti*PJ300WY01	PT+PU400WY
AFF clip	Yes	Yes	Yes
Filter (µm)	15	15	15
Pump tube (mm)	2000	2000	2000
Injector port	None	Needle port	Needle port
Drip rate (drops/ml)	20	20	20
Material - tubing	PVC ¹	Polybutadiene ⁴	PVC ¹
Packaging (unit/carton)	20/200	20/200	25/250
Priming volume (ml) ²	15.7	16.7	18.1
Colour ³	Transparent	Transparent	Transparent
Pressure resistance (bar)	1.5	1.5	1.5
Total length of sets	2150	2440	2380
Product Illustration			

General (amber)

Ti*PA300WY01E	Ti*PA300WN01E
Yes	No
15	15
2000	2000
Needle port	Needle port
20	20
PVC ¹	PVC ¹
20/200	20/200
18.2	18.2
Amber	Amber
1.5	1.5
2640	2640

Transfusion

TB*PU300SYE
Yes
200
2200
None
20
PVC ¹
20/200
17.2
Transparent
1.5
2420

Specific

SA*QTF301NN	SA*QNT310MN	SA*QNT310YN	SA*QNT320MN
Yes	No	Yes	No
0.2	40	40	40
1220	1220	1220	1220
T-Port ⁵ x 1	3WSC ⁶ x 1	3WSC ⁶ x 1	3WSC ⁶ x 2
20	20	20	20
Polybutadiene ⁴	Polybutadiene ⁴	Polybutadiene ⁴	Polybutadiene ⁴
20/160	20/160	20/160	20/160
15.4	13.5	11.8	13.9
Transparent	Transparent	Transparent	Transparent
1.5	1.5	1.5	1.5
2450	2400	1900	2450

Specific

Product code	SA*QNT320YN	SA*QNT330MN	SA*QNT330YN
AFF clip	Yes	No	Yes
Filter (µm)	40	40	40
Pump tube (mm)	1220	1220	1220
Injector port	3WSC ⁶ x 2	3WSC ⁶ x 3	3WSC ⁶ x 3
Drip rate (drops/ml)	20	20	20
Material - tubing	Polybutadiene ⁴	Polybutadiene ⁴	Polybutadiene ⁴
Packaging (unit/carton)	20/160	20/160	20/160
Priming volume (ml) ²	12.2	14.3	14.3
Colour ³	Transparent	Transparent	Transparent
Pressure resistance (bar)	1.5	1.5	1.5
Total length of sets	1950	2500	2500
Product Illustration			

SA*QTF301UN	SA*QTF302NN	SA*QTF302UN	SA*QTF303NN	SA*QTF303UN
No	Yes	No	Yes	No
0.2	0.2	0.2	0.2	0.2
1220	1220	1220	1220	1220
T-Port ⁵ x 1	T-Port ⁵ x 2	T-Port ⁵ x 2	T-Port ⁵ x 3	T-Port ⁵ x 3
20	20	20	20	20
Polybutadiene ⁴	Polybutadiene ⁴	Polybutadiene ⁴	Polybutadiene ⁴	Polybutadiene ⁴
20/160	20/160	20/160	20/160	20/160
15.4	15.1	15.1	15.8	15.8
Transparent	Transparent	Transparent	Transparent	Transparent
1.5	1.5	1.5	1.5	1.5
2450	2350	2350	2500	2500

SA*QTT301LN	SA*QTT301MN	SA*QTT301YN	SA*QTT302MN	SA*QTT302YN
No	No	Yes	No	Yes
15	40	40	40	40
1220	1220	1220	1220	1220
T-Port ⁵ x 1	T-Port ⁵ x 1	T-Port ⁵ x 1	T-Port ⁵ x 2	T-Port ⁵ x 2
20	20	20	20	20
Polybutadiene ⁴	Polybutadiene ⁴	Polybutadiene ⁴	Polybutadiene ⁴	Polybutadiene ⁴
20/160	20/160	20/160	20/160	20/160
11.7	11.7	11.7	12.2	12.2
Transparent	Transparent	Transparent	Transparent	Transparent
1.5	1.5	1.5	1.5	1.5
1900	1900	1900	2050	2050

1 Does not contain DEHP - 2 Total volume of the infusion line (drip chamber empty) - 3 Tubing & Drip chamber - 4 Does not contain PVC - 5 SURPLUG AD T-Port Needleless connecting device - 6 SURPLUG AD 3-way stopcock

Terufusion - Full focus infusion system



Terumo offers a comprehensive portfolio of solutions for the complete infusion line-up. User-friendly syringe and volumetric pumps and racks, as well as a wide range of highly accurate, easy-to-use disposables, such as an extended range of infusion sets and syringes.



TERUMO EUROPE N.V.

Terumo Corporation Terumo Europe N.V.
+81 3 6742 8500 +32 16 38 12 11

www.terumo-europe.com

Your Terumo partner

This brochure is only applicable for countries within EMEA, excluding the MENA region.

Please note that some products may not be available in all countries, check with your local Terumo contact.

All brand names are trademarks or registered trademarks owned by TERUMO CORPORATION, its affiliates, or unrelated third parties. Published by Terumo Europe N.V.