

DF55

Blockbuster Mini 5-Part Hematology Analyzer

Principles

Flow cytometry(FCM)+Tri-anglescatter for 5-part differential analysis and WBC counting
Impedance method for RBC and PLT counting
Cyanide free colorimetry for HGB test

Parameters

25 reportable parameters: WBC, Neu#, Lym#, Mon#, Eos#, Bas#, Neu%, Lym%, Mon%, Eos%, Bas%, RBC, HGB, HCT, MCV, MCH, MCHC, RDW-SD, RDW-CV, PLT, MPV, PDW, PCT, P-LCR, P-LCC
6 research parameters: ALY#, ALY%, LIC#, LIC%, NRBC#,NRBC%
3 histograms for WBC, RBC and PLT
Three DIFF scattergrams and one BASO scattergram for WBC differentiation

Sample volume

20 μ L

Linearity range

WBC: 0.00-300.00 $\times 10^9$ /L
RBC: 0.00-8.50 $\times 10^{12}$ /L
HGB: 0-250g/L
PLT: 0-3000 $\times 10^9$ /L
HCT: 0.0-67.0%

Repeatability

WBC $\leq 2.0\%$ (4.0-15.0 $\times 10^9$ /L)
RBC $\leq 1.5\%$ (3.5-6.0 $\times 10^{12}$ /L)
HGB $\leq 1.5\%$ (110-180g/L)
MCV $\leq 1\%$ (70-120fL)
PLT $\leq 6.0\%$ (100-149 $\times 10^9$ /L)
 $\leq 4.0\%$ (150-500 $\times 10^9$ /L)

Sample mode

Whole blood, capillary whole blood and pre-diluted modes

Throughput

Up to 60 tests per hour

Display

10.4 inches TFT touch screen

Printout

Internal printer: 57 mm thermal printer with various formats
External printer: compatible with multiply laser / inkjet printers, compatible with various formats and user-defined formats

Storage

Up to 100,000 records

Communication

LAN port supports HL7 protocol
Support bi-directional LIS

Interface

4 USB ports, 1 LAN port

Power requirement

100V-240V~, 50/60Hz, ≤ 200 VA


Dimension

364mm(W)*498mm(H)*431mm(D)

Net weight

28kg

Expert for Small Lab



10th Floor, Building B, High-tech Park, Guangqiao Road, Tianliao Community, Yutang Street, Guangming District, Shenzhen 518107, P. R. China

+86-755-26008015
www.dymind.com

+86-755-26008015
intl@dymind.com



Declaration: Shenzhen Dymind Biotechnology Co., Ltd reserves the right to change the product of specifications and appearance at any time. For the information of this manual, Shenzhen Dymind Biotechnology Co., Ltd reserves the right to the interpretation and the decision.
P/N: EN-DF55[5.0]

DF55

Blockbuster Mini 5-Part Hematology Analyzer

DF55

Blockbuster Mini 5-Part Hematology Analyzer

Tri-anglelaser scatter with advanced flow cytometry method

Tri-angle semi-conductor laser scatter can show cellular complexity, cellular granularity and cell volume, that means we can get accurate WBC differentiation.

Advanced flow cytometry can get stable sheath fluid and ensure every cell goes through the flow cell one by one.

Independent BASO channel

Basophils (BASO) has important clinical significance, such as Leukemia, Anaphylactic Disease, Hematemesis, Cancer and so on.

Real double optical channel test, both for DIFF and BASO, independent BASO channel with optical counting contributes to more precise results.

Easy to operate

10.4 inches TFT touch screen with 12.5° best angle

Saving space

Built-in reagent room is convenient for reagent positioning

Compact size

Small footprint for space saving

Rich interface

4 USB ports, 1 LAN port

Humanized design

No barrier around the needle,
It is suitable for either left hand or right hand

