

PROVIEW 10 / 12 Accessory list





Table of Contents

• Standard Accessory 2	
• ECG	3
• SpO2	
• Non-Inva	sive Blood Pressure 5
• Temperat	ure
• Invasive B	Blood Pressure 6
• C.O	
• ETCO2	8
MicroFlo	w Stream (BLT)
Main Str	eam (BLT)
• Others	

1

Standard accessory



5-lead Snap Connector ECG Cable

41.20-5362

IEC, 12 pin, AQS, Encrzpted
Compatible with PROVIEW10, 12



Adult Disposable Adhesive Electrodes

70.10-3700

Snap Connector, 30pcs / pack
Compatible with PROVIEW10, 12



Analog SpO2 adult soft reusable sensor

41.20-5363

12 pin, Encrypted, One piece Compatible with PROVIEW10, 12



NIBP PVC hose

41.20-5369

3 m, A series connector

Compatible with PROVIEW10, 12



Adult cuff

21.10-7051

with connector M5125, 250mm \sim 350 mm Compatible with PROVIEW10, 12



Lithium Battery

41.20-5375

11.1 V / 2500mAh

Compatible with PROVIEW10, 12



User's Manual

Manual

220 V

Compatible with PROVIEW10, 12

ECG



5-lead Snap Connector ECG Cable

41.20-5362

IEC, 12pin, AQS, Encrypted
Compatible with PROVIEW10, 12



3-lead Snap Connector ECG Cable

41.20-5361

IEC, 12pin, AQS, Encrypted
Compatible with PROVIEW10, 12



Adult Disposable Adhesive Electrodes

70.10-3700

Snap Connector, 30 pcs/pack Compatible with PROVIEW10, 12




Pediatric / Neonatal Disposable Adhesive Electrodes


70.10-2107

Snap Connector, 30 pcs/pack
Compatible with PROVIEW10, 12

SpO₂



NIBP



TEMP



Skin Temperature Probe

for adult, pediatric, neonatal
Compatible with PROVIEW10, 12

41.20-5371



Rectal Temperature Probe

for adult, pediatric, neonatal

Compatible with PROVIEW10, 12

41.20-5372

IBP



IBP Kit

Disposable pressure transducer kit (IBP cable 1pc, IBP disposable transducer 25 pcs) Compatible with PROVIEW 12 41.20-5389



IBP Cable

with connector

Compatible with PROVIEW 12

41.20-5388



IBP Disposable Transducer

Compatible with PROVIEW 12

21.11-5951





C.O. Module

41.20-5381

without Accessory Kit

Compatible with PROVIEW 12



C.O. Kits

41.20-5382

- 1 Cardiac Output cable
- In-line Injection temperature probe (SP4042)
 In-line Injection temperature probe housing (SP5045)
- 1 Control Syringe

Compatible with PROVIEW 12



Cardiac Output cable

41.20-5441

with connector

Compatible with PROVIEW 12



In-line Injection temperature probe

41.20-5442

SP4042

Compatible with PROVIEW 12



In-line Injection temperature probe housing

41.20-5443

BD SP5045

Compatible with PROVIEW 12



Control Syringe

41.20-5444

Compatible with PROVIEW10, 12

ETCO2

MicroFlow Stream



Capno_S ETCO2

41.20-5385

Water filter 2 pcs, Disposable CO2 sample tube 1pc, Airway adapter 1pc Compatible with PROVIEW 12



12 Pin ETCO2 Module Extension Cable

41.20-5390

for Masimo Gas Module

Compatible with PROVIEW 12



Water Filter

41.20-5445

T4F

Compatible with PROVIEW 12



L type 3 way airway adapter

41.20-5446

for intubated patient

Compatible with PROVIEW 12



Disposable CO2 Nasal Cannula

41.20-5447

for Adul

Compatible with PROVIEW 12

ETCO2

Main Stream



Capno_M ETCO2

41.20-5386

with Disposable airway adapter for Adult / Pediatric 1pc Compatible with PROVIEW 12



12Pin ETCO2 module extension cable

41.20-5390

for Masimo Gas module

Compatible with PROVIEW 12




Disposable Airway Adapter

41.20-5448

for Adult / Pediatric

Compatible with PROVIEW 12

Others



ver. 002

Distributor

*All specifications are subject to change without notice

Copyright @ 2020 medical ECONET All rights reserved