


DAHLHAUSEN®
Lotus 15


Laryngeal Mask, single use

The sterile Laryngeal Mask, single-use, provides a respiration without intubation of the trachea. It seals the larynx to ensure safe respiration.

Material: Silicone, Polypropylene, Polycarbonat

Sterile: Yes

Art.no.	Application area	Max. filling volume (cuff)	Maximum body weight	Minimum body weight
4400060010	Newborn	4 ml	5 kg	<5 kg
4400060015	Small Child	7 ml	10 kg	5 kg
4400060020	Child	10 ml	20 kg	10 kg
4400060025	Child	14 ml	30 kg	21 kg
4400060030	Adolescent	20 ml	50 kg	31 kg
<u>4400060040</u>	<u>Adult</u>	<u>30 ml</u>	<u>70 kg</u>	<u>51 kg</u>
<u>4400060050</u>	<u>Adult</u>	<u>40 ml</u>	<u>100 kg</u>	<u>71 kg</u>

P. J. Dahlhausen & Co. GmbH
Emil-Hoffmann-Str. 53
50996 Köln