Product Specifications

Tissue-Tek[®] Mounting Media





Tissue-Tek[®] Glas[™] Pertex

Product description

Item Code	Product Description	Sales quantity
1410	Glas™ Pertex	6x500 mL

Product specifications

1. Application	
	Mounting medium for coverslipping glass slides

2. Related Instrument		
Item Code	Product Description	
6502	Glas [™] g2 Glass Coverslipper	
3. Regulatory information		
Intended Use	Mounting media work as a glue to attach coverslips to glass slides to	

		factors, while being as see-through as the glass itself.	
CE status			
MSDS available	Yes ⊠(14 Languages)	No 🗆	
CoA available	Yes 🛛	No 🗆	

Product Specifications Tissue-Tek[®] Mounting Media

4. Packaging		
Package description	Box 6 bottles	
Sales quantity	1410	6 of 500 mL
Dimensions (LxWxH) or (ØxH)	76x155mm	
Package dimensions (LxWxH)	268x178x186mm	
Shipping dimensions (LxWxH)	268x178x186mm	
Shipping weight	4 kg	

5. Labeling		
HIBC code	+H905140821	
EAN code	NA	
Statistical code	32082010	



6. Storage and transport		
Storage conditions:		
 Temperature 	18 – 30°C	
Humidity	Dry environment	
 Sunlight 	Protect from sunlight	
Transport conditions:		
 Temperature 	Frost-free and max. at 35°C	
 Humidity 	Dry environment	
 Sunlight 		
Handle with care (fragile)	Yes 🛛	No 🗆
Dangerous Good	Yes 🛛	No 🗆
Shelf life	Expiry date on product	

7. Material specifications		
Components	Xylene (25 – 50%), Ethylbenzene (10 – 20%)	
Chemical resistance	Keep away from oxidizing agents	
Heat resistance	Store frost free in a cool place, protect from sunlight	
Cleaning	Collect spilled materials in container, absorb residues with sand or other inert material.	
Latex free	Yes* 🗆 No 🛛	
	*only when official statement is available	

8. Terms & Conditions

Product Specifications Tissue-Tek[®] Mounting Media

Warranty	6 months after receiving by buyer
Terms of Delivery	F.O.B. Warehouse, The Netherlands
9. Additional information	