

Blood Collection Tubes Specification Sheet

Category	Colour	Additive	Specification	Volume	Material	Package
Serum Blood Collection Tubes		No	13×75mm	1~5ml	Glass/PET/PP	100PCS/TRAY 1800PCS/CTN
			13×100mm	5~7ml		
		Clot Activator	13×75mm	1~5ml	Glass/PET/PP	
			13×100mm	5~7ml		
		Gel&Clot Activator	13×75mm	1~5ml	Glass/PET/PP	
			13×100mm	5~7ml		
Plasma Blood Collection Tubes		Sodium Fluoride	13×75mm	1~5ml	Glass/PET/PP	
			13×100mm	5~7ml		
		Sodium Citrate 9NC	13×75mm	1~5ml	Glass/PET/PP	
			13×100mm	5~7ml		
		Sodium Heparin Lithium Heparin	13×75mm	1~5ml	Glass/PET/PP	
			13×100mm	5~7ml		
Blood Collection Tubes		EDTA.K2 EDTA.K3	13×75mm	1~5ml	Glass/PET/PP	
			13×100mm	5~7ml		
		Sodium Citrate 4NC	13×75mm	1~5ml	PET/PP	
			13×100mm	5~7ml		
			8×120mm	1.28ml	Glass	100PCS/TRAY 1200PCS/CTN
			8×120mm	1.6ml		

Mini Blood Collection Tubes Specification Sheet

Category	Colour	Additive	Specification	Volume	Material	Package
Serum Blood Collection Tubes		No	8x44mm	0.2~1ml	PP	100PCS/TRAY 4800PCS/CTN
		Clot Activator	8x44mm	0.2~1ml		
		Gel&Clot Activator	8x44mm	0.2~1ml		
Plasma Blood Collection Tubes		Sodium Fluoride	8x44mm	0.2~1ml		
		Sodium Heparin Lithium Heparin	8x44mm	0.2~1ml		
Blood Collection Tubes		EDTA.K2 EDTA.K3	8x44mm	0.2~1ml		

Blood Collection Needles

Butterfly type

It is mainly used with Vacuum Blood Collection Tubes for blood specimen collection in the clinical laboratory.

REF	Specification	Color
TYPE-A-20G	□ 20G	Yellow
TYPE-A-21G	□ 21G	Green
TYPE-A-22G	□ 22G	Black
TYPE-A-23G	□ 23G	Blue
TYPE-A-25G	□ 25G	Orange



Butterfly type with safety device

It is mainly used with Vacuum Blood Collection Tubes for blood specimen collection in the clinical laboratory.

REF	Specification	Color
TYPE-B-21G	□ 21G	Green
TYPE-B-22G	□ 22G	Black
TYPE-B-23G	□ 23G	Blue
TYPE-B-25G	□ 25G	Orange



Pen type

It is mainly used with Vacuum Blood Collection Tubes for blood specimen collection in the clinical laboratory.

REF	Needle Specification		Color
TYPE-D-18G	□ 18GX1"	□ 18GX1 1/2"	Pink
TYPE-D-20G	□ 20GX1"	□ 20GX1 1/2"	Yellow
TYPE-D-21G	□ 21GX1"	□ 21GX1 1/2"	Green
TYPE-D-22G	□ 22GX1"	□ 22GX1 1/2"	Black
TYPE-D-23G	□ 23GX1"	□ 23GX1 1/2"	Blue



Pen type with safety device

It is mainly used with Vacuum Blood Collection Tubes for blood specimen collection in the clinical laboratory.

REF	Needle Specification		Color
TYPE-C-18G	□ 18GX1"	□ 18GX1 1/2"	Pink
TYPE-C-20G	□ 20GX1"	□ 20GX1 1/2"	Yellow
TYPE-C-21G	□ 21GX1"	□ 21GX1 1/2"	Green
TYPE-C-22G	□ 22GX1"	□ 22GX1 1/2"	Black
TYPE-C-23G	□ 23GX1"	□ 23GX1 1/2"	Blue

