










































| Compatible Brand | P/N     | Description  | Male Connector  | Skin/Oral Side  |
|------------------|---------|--|---|---|
| Comen            | U638-1D | Pediatric Skin Temperature probe,10KΩ,2pin,3meter                      |    |    |
|                  | U638-1N | Infant Skin Temperature probe,10KΩ,2pin,3meter                         |    |    |
| Creative         | U610-1A | Adult Rectal/Esophageal Temperature probe,10KΩ,2pin,3meter             |    |    |
|                  | U610-1C | Pediatric Rectal/Esophageal Temperature probe,10KΩ,2pin,3meter         |   |    |
|                  | U610-1M | Infant Rectal/Esophageal Temperature probe,10KΩ,2pin,3meter            |   |    |
|                  | U610-1B | Adult Skin surface Temperature probe,10KΩ,2pin,3meter                  |   |   |
|                  | U610-1D | Pediatric Skin Temperature probe,10KΩ,2pin,3meter                      |   |  |
|                  | U610-1N | Infant Skin Temperature probe,10KΩ,2pin,3meter                         |   |  |
| Datascope        | U614-1A | Datascope Adult Rectal/Esophageal Temperature probe,2.25KΩ,5pin,3meter |  |  |
|                  | U614-1C | Rectal/Esophageal Temperature probe,2.25KΩ,5pin,3meter                 |   |  |
|                  | U614-1M | Infant Rectal/Esophageal Temperature probe,2.25KΩ,5pin,3meter          |  |  |
|                  | U614-1B | Adult Skin surface Temperature probe,2.25KΩ,5pin,3meter                |   |  |
|                  | U614-1D | Pediatric Skin Temperature probe,2.25KΩ,5pin,3meter                    |   |  |

| Compatible Brand         | P/N  | Description   | Male Connector  | Skin/Oral Side  |
|--------------------------|--|---|---|---|
| Datascope                | U614-1N  | Infant Skin Temperature probe,2.25KΩ,5pin,3meter                                      |    |    |
| Drager/Air Shields C2000 | U615-2N  | Drager Infant Skin Temperature Probe,Mini5Pin,9.1KΩ,3meter                            |    |    |
| GE                       | U602-1A  | Adult Rectal/Esophageal Temperature probe,2.25KΩ,11pin,3meter                         |    |    |
|                          | U602-1C  | Pediatric Rectal/Esophageal Temperature probe,2.25KΩ,11pin,3meter                     |   |    |
|                          | U602-1M  | Infant Rectal/Esophageal Temperature probe,2.25KΩ,11pin,3meter                        |   |    |
|                          | U602-1B  | Adult Skin surface Temperature probe,2.25KΩ,11pin,3meter                              |   |   |
|                          | U602-1D  | Pediatric Skin Temperature probe,2.25KΩ,11pin,3meter                                  |   |  |
|                          | U602-1N  | Infant Skin Temperature probe,2.25KΩ,11pin,3meter                                     |   |  |
|                          | U602-2C  | GE V100 Oral Temperature Probe,Crystal head Connector                                 |  |  |
|                          | U602-1F  | GE Temperature interface cable,Rectangle 2Pin,3meter                                  |  |  |
|                          | U602-1G  | GE Temperature Interface Cable,3.5mm Female   |  |  |
|                          | U602-1J  | GE Temperature extension cable,to YSI 400 series & YSI 700 series                     |  |  |
| U602-3J                  | GE Temperature interface cable,6.3mm Female,0.5meter |  |  |   |