



**GT029-800**

**Disposable Colostomy Bag**

Size: 20cm x 14cm, 25cm x 15cm  
 Adhesive area: 10cm x 9cm, 12cm x 12cm  
 Center opening:  $\phi$ 3cm  
 Sterilization: EO  
 Form of packing:  
 20pc/PE bag, 5 bags/box, 20boxes/carton  
 Measurement of carton:  
 470x330x240mm(20cm x 14cm)  
 580x350x240mm(25cm x 15cm)



**GT029-900V**

**Disposable Enema Bag**

Capacity: 1600ml  
 Material: PVC  
 Sterilization: EO  
 Form of packing:  
 1pc/PE bag, 200pcs/carton  
 Measurement of carton:  
 460x440x330mm



**GT031-100**

**Disposable Enteral Feeding Bag Set**

Size: 500ml, 1000ml, 1200ml  
 Available for Pump Using and Gravity Using  
 Packed in individual PE Bag  
 30pcs/ctn,  
 Dimension: 40X28X25cm

**GT031-500**

**Irrigation syringe**

Fiat: 80ml, 60ml  
 Parts: Protector cap(PE), barrel(PP), piston(synthetic rubber), handspike(PP), two way(PP) 400pcs/58x44x46cm Ring: 60ml  
 Parts: Protector cap(PE), barrel(PP), piston(synthetic rubber), handspike(PP), two way(PP) 250pcs/58x44x46cm  
 Bulb: 60ml, 60ml Parts: Protector cap(PE), barrel(PP), bulb(PE), two way(PP) 50pcs/35x26x34cm



**NEW STYLE OSTOMY BAG**

Our ostomy bag are of excellent quality and advanced usage, making the wound healed in the good condition. Ostomy bag applications: collect excrement and urine after colostomy, ileostomy and urostomy.

**Features:**

1. The hydrocolloid skin barrier has the performance of anti-maceration, low allergy, and warping and leakage. It has excellent exudation—absorbent and long sticking ability(up to 7days). And it is soft and comfortable.
2. Hydrocolloid skin barrier combines with the highly permeable PU film makes it stick for a long time, no leakage of exudates, convenient to use.
3. The material of ostomy bag possesses high obstructing ability, it is comfortable and covert.
4. The one—piece system design discharge intestinal gas easily and filtering effect it excellent.

Simple type  
 Foam type  
 Glue Pan type  
 Two—piece type

Packing Reference:  
 10pca/bag, 400pca/ctn  
 Size: 350x310x400mm  
 G.W: 6kg N.W: 4.5kg



**GT030-200 Simple Type**



**GT030-300 Foam Type**



**GT030-400 Glue Pan Type**



**GT030-500 Two—piece Type**



**GT030-600 Foam Type**



**GT030-700 Two—piece Type**