

# Mollelast® haft

## Cohesive Bi-Elastic Conforming Bandage

### One... and done



The cohesive conforming bandage with **two-way stretch** and special expansion behaviour adheres to itself, but not to skin or hair. It is available in different versions for a variety of needs: Mollelast haft, Mollelast haft color and Mollelast haft latex-free.

Mollelast haft stretches lengthwise and widthwise, so there is no creasing and it stays in place easily. Can be applied easily to tapered parts of the body and even stays in place reliably around joints thanks to the excellent adhesion between layers. This product can also be counted on to stay in place during movement: no constriction or slipping.

The good cohesive properties come from a microcoating. The selvages make application a breeze and prevent fraying. The thin knitted material is very permeable to air and gentle to the skin.

The bandage is used for the fixation of dressing pads, cannulas and splints (such as Cellacast Splint) and can provide support, for example in the case of minor sports injuries. Can be used in combination with Cellacast Splint.

### Latex-free fixation

Mollelast haft also is available latex-free, i.e. with a synthetic microcoating (free from natural rubber latex). With the latex-free product, users have as an option a very well-tolerated conforming bandage which has all the tried-and-tested properties of Mollelast haft.

### Product composition

**Mollelast haft, 66 % viscose, 34 % polyamide, coated with microfine natural rubber adhesive (latex).** Mollelast haft color, 56 % viscose, 44 % polyamide, coated with microfine natural rubber adhesive (latex)

## Versions

Mollelast haft: lengthwise elastic to about 65%/widthwise elastic to about 50%

Mollelast haft color: lengthwise elastic to about 55%/widthwise elastic to about 50%

## Ordering information

Mollelast haft color

Mollelast haft

stretched length, individually packed in a folding box

Dimensions (cm x m)	REF	SC/TC
<b>white</b>		
4 x 4	30063	1/120
6 x 4	30064	1/100
8 x 4	30065	1/120
10 x 4	30066	1/120
12 x 4	30067	1/100
4 x 20	30068	1/18
<b>Lotul 2</b> 6 x 20	30069	1/18
8 x 20	30070	1/18
<b>Lotul 1</b> 10 x 20	30071	1/18
12 x 20	30072	1/18

## Fields of application

- for immobilisation of limbs
- for retention of dressings and cannulas

## Properties

- lengthwise elastic
- widthwise elastic
- cohesive because of latex particles
- air-permeable
- is not very bulky
- the turns of the bandage stay reliably in place
- economical use of materials
- high-quality, stable selvages