

Transfer Bag System




- **SAFE BLOOD**

- **SAFE PRODUCT**

Mitra Blood Bags are designed for the collection, processing and storage of whole human blood and it's components.

...complete solution for **BLOOD** transfusion services

Transfer Bags are used for preservation and transfusion of Whole Human Blood or it's Components.



1. Spike Cover
2. Spike assembly
3. PVC Tube
4. External Port
5. Side Slits
6. Label
7. Bottom Slit

TECHNICAL SPECIFICATIONS

- Transfer Bag Plain - 100 / 150 / 200 / 250 / 300 / 500 / 600 / 800 / 1000 / 2000 ml.
- DEHP plasticized PVC, collapsible, non-vented, sterile container complete with collecting tube in a closed system.
- Sterile non-pyrogenic, non-toxic, non-hemolytic, bio-compatible material.
- Inner bag maintains the sterility.
- No risk of contamination & air embolism with all leak proof.
- Slits on the both side to accommodate 5-10 ml volume test tube & at the bottom of bag to hang the blood bag during transfusion.
- Flexible, non-sticking, transparent, leak proof tubing having multiple printed & ID/Segment number.
- Polycarbonate sharp spike, rust proof, tightly fixed with tube.
- Spike assembly & easily accessible external ports.
- Hot melt technology non-peel able, heat sealed label with transparent adhesive remarked attached between room temperature to -80°C.
- 3 years shelf life from the date of manufacturing.
- Protective dual packaging (CPP & Corrugated Box)
- Blood bags meet all the standards as laid down in ISO 3826.
- Individual plastic transfer bags packed in plastic CPP pouch & such 25 plastic pouches packed in duplex box & 4 duplex box packed in corrugated box.

Specification/Designs are subjected to change without any obligation on the part of the manufacturer

WHO GMP
CERTIFIED

CE
2460



Manufactured & Marketed by:

Mitra Industries Pvt. Ltd.

14/4, Delhi-Mathura Road, Faridabad, Haryana-121 003. INDIA
T: +91-129-2271068/69, 4099300 | Toll Free: 1800 1020 307
E: customersupport@mipl.in | W: www.mitraindustries.com

MIPL/11/12

Sistem de sac de transfer

Siguranta sangvina .

Produs in siguranta.

Pungile de sânge Mitra sunt concepute pentru colectarea, prelucrarea și depozitarea sângelui uman întreg și sunt componente.

soluție completă pentru transfuzie și vicii de sânge...

Pungile de transfer sunt folosite pentru conservare și transfuzie de sânge uman întreg sau componentele sale.

1. Capacul Spike
2. Ansamblu spike
3. Tub din PVC
4. Port extern
5. Fante laterale
6. Eticheta
7. Fanta de jos

SPECIFICATII TEHNICE

Pungă de transfer simplă - 100/150/200/250/300/500/600/800/1000/2000 ml.

- PVC plastificat DEHP, container pliabil, fără aerisire, steril, complet cu tub de colectare într-un sistem închis.
- Material steril non-pirogen, netoxic, non-hemolitic, bio-compatibil.
- Geanta interioară menține sterilitatea.
- Niciun risc de contaminare și embolie aeriană cu toate etanșeitățile.
- Fante pe ambele părți pentru a găzdui o probă de volum de 5-10 ml și în partea inferioară a pungii pentru a agăța punga de sânge în timpul transfuziei.
- Tuburi flexibile, antiaderente, transparente, rezistente la scurgeri, cu mai multe numere tipărite și de identificare / segment.
- Vârf ascuțit din policarbonat, rezistent la rugină, bine fixat cu tub.
- Asamblare vârf și porturi externe ușor accesibile
- Tehnologie de topire la cald, etichetă sigilată termic, care nu se poate dezlipi, cu adeziv transparent
- 3 ani de valabilitate de la data fabricației.
- Ambalaj dublu de protecție (CPP și cutie ondulată).
- Pungile de sânge îndeplinesc toate standardele prevăzute în ISO 3826.
- Pungi individuale de transfer din plastic ambalate în pungă din plastic CPP și astfel de 25 de pungi din plastic ambalate în cutie duplex și 4 cutii duplex ambalate în cutie ondulată.