

## SANAVEL

### Dressing fabric SANAVEL

Hydrophobic high-porous dressing fabric SANAVEL is a warp-knitted fabric made of polyester filament yarn. SANAVEL 1 differs in its knitted structure from SANAVEL 2.

#### Indications:

The SANAVEL dressing is used for treatment of both open and closed wounds. Impregnation with drugs in a liquid or ointment base prior to application is advisable. It is designed particularly for use in surgical and dermatological treatment, in clinical plastic surgery, in both hospitals and ambulatory care facilities.

SANAVEL 1 - is suitable for impregnation with drugs in an ointment base

SANAVEL 2 - is suitable for impregnation in a liquid vehicle



- it is non-adhesive, has a pleasant feel, does not irritate the wounds and does not stick to them
- it can be easily removed from any wound without the risk of damage to the wound base
- it displays good secretion permeability, which favorably influences the healing processes
- the expiration period is 3 years

## Product types, sizes and catalogue codes:

### Dressing fabric SANAVEL

Type	Size [mm]	Number of pieces in package	Catalogue number
SANAVEL 1	140 × 100	25 pcs	131-1410
SANAVEL 2	140 × 100	25 pcs	132-1410

*Other dimensions are produced on the basis of customer specifications.*

### **VUP Medical, a.s.**

Šujanovo nám. 3 | 602 00 Brno

Czech Republic | [www.vup.cz](http://www.vup.cz)

phone: +420 543 160 444

e-mail: [sales@vup.cz](mailto:sales@vup.cz)