

TECHNICAL DATA SHEET

AUTOCLAVABLE BIOHAZARD BAG

PRODUCT CODE	:	OT
RELATED STANDARDS	:	Autoclavable Biohazard Bag meets the requirements of the Council Directive 94/62/EC on Packaging.
INTENDED USE	:	Autoclavable Biohazard Bags are used for the sterilization of medical wastes before disposal.
OPERATION CONDITIONS	:	Autoclavable Biohazard Bags are designed for use in steam sterilization methods. The recommended operation temperature is max. 141 °C. Other sterilization conditions should be determined by the end user regarding to material to be sterilized. During sterilization, bag should be left open for an effective sterilization.

PRODUCT SPECIFICATIONS

PROPERTIES	METHOD	UNIT	VALUE	
THICKNESS (one side)	ISO 534	µm	50 ± 5	
TENSILE STRENGTH	ISO 527-3	N/15 mm	MD	min. 12
			CD	min. 11
ELONGATION AT BREAK	ISO 527-3	%	MD	min. 50
			CD	min. 30
TEAR FORCE	ASTM D 1922	mN	min. 700	

PACKAGING : Autoclavable Biohazard Bags are packaged as the listed below.

<u>Product Code</u>	<u>Dimensions</u>	<u>Unit</u>
OT5370	53 cm x 70 cm	200 pcs/box
OT6389	63 cm x 89 cm	120 pcs/box
OT8790	87 cm x 90 cm	100 pcs/box
OT87110	87 cm x 110 cm	100 pcs/box
OT80110	80 cm x 110 cm	100 pcs/box



PMS Tıp Teknolojileri San. ve Tic.A.Ş
pms@pmsmedikal.com