

CHIRALEN monofilament

UK07
142



Non-absorbable - Chiralen monofilament

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
2 x HR 8 	0.5	7/0	1 x 45 cm	Blue	PP 5205	
HR 10 	0.7	6/0	1 x 45 cm	Blue	PP 5142	
2 x HR 10 	0.7	6/0	1 x 45 cm	Blue	PP 5210	
	1	5/0	1 x 45 cm	Blue	PP 5211	
	1	5/0	1 x 60 cm	Blue	PP 6211	
HR 13 	0.7	6/0	1 x 75 cm	Blue	PP 5145	
	1	5/0	1 x 75 cm	Blue	PP 5146	
2 x HR 13 	1	5/0	1 x 45 cm	Blue	PP 5220	
	0.7	6/0	1 x 75 cm	Blue	PP 5219	
	1	5/0	1 x 75 cm	Blue	PP 5222	
2 x HR 15 	1.5	4/0	1 x 90 cm	Blue	PP 5228	
	1	5/0	1 x 90 cm	Blue	PP 5558	
HR 17 	1	5/0	1 x 90 cm	Blue	PP 5558	
2 x HR 17 	1	5/0	1 x 75 cm	Blue	PP 5236	
	1	5/0	1 x 90 cm	Blue	PP 5248	
HR 18 	1.5	4/0	1 x 75 cm	Blue	PP 5153	
	2	3/0	1 x 75 cm	Blue	PP 5154	
2 x HR 18 	1	5/0	1 x 75 cm	Blue	PP 5233	
	1.5	4/0	1 x 75 cm	Blue	PP 5234	
	2	3/0	1 x 75 cm	Blue	PP 5235	
2 x HR 18 	1	5/0	1 x 90 cm	Blue	PP 5237	
	1.5	4/0	1 x 90 cm	Blue	PP 5238	
	2	3/0	1 x 90 cm	Blue	PP 5239	

CHIRALEN monofilament

UK07
141

CHIRALEN monofilament is a non-absorbable suture made of polypropylene. The thread is hydrophobic and is not subject to degradation by tissue enzymes.

Chiralen monofilament has a smooth, non-porous surface, that guarantees a clean and smooth passage through the tissue, optimal elasticity, good elongation properties and high visibility. In tissues, it can cause a minimal acute inflammatory reaction.



Indication: All surgical specialities, especially Cardiovascular, Plastic, Reconstructive, Ophthalmic, Orthopaedic and Traumatology (e.g. tendons).

Contraindication: Chiralen monofilament must not penetrate into the inner parts of the digestive and urinary tracts.

Chemical Composition: Polypropylene

Sizes: 10/0 to 2 USP

Structure: Monofilament

Absorption method: Non-absorbable

Colour: Blue and Black

Non-absorbable - Chiralen monofilament