


AlCath G FullCircle

4 mm GoldTip Ablation Catheter



Ordering Information

Model	Tip	Curve	Spacing	Shaft Diameter	Length	Order Number
AlCath Red TC G FullCircle	4 mm	Red	2-5-2 mm	7 F	110 cm	364537
AlCath Blue TC G FullCircle	4 mm	Blue	2-5-2 mm	7 F	110 cm	364525
AlCath Green TC G FullCircle	4 mm	Green	2-5-2 mm	7 F	110 cm	364533
AlCath Black TC G FullCircle	4 mm	Black	2-5-2 mm	7 F	110 cm	364522
AlCath Cyan TC G FullCircle	4 mm	Cyan	2-5-2 mm	7 F	110 cm	364529

Product Highlights

GoldTip Ablation

GoldTip Ablation results in 33% faster lesions compared to conventional catheters¹

FullCircle Deflection

Deflectability beyond 180 degrees through FullCircle deflection²

Handling

Precise positioning through push-pull steering incl. AutoLock²

¹ T. Balazs et al., J Interv Card Electrophysiol. 2012 Oct 26; DOI: 10.1007/s10840-012-9743-9.

² VAL11_AlCath_019_Dimensional Test_PG0_Rev-0

VAL11_AlCath_019_Dimensional Test_PG4_Rev-0

AlCath G FullCircle

Technical Data

Catheter Specification

Tip electrode length	4 mm
Tip electrode type	Non-irrigated
Tip temperature sensor	Thermocouple
Ring electrode length	1.55 mm
Number of electrodes	4
Electrode spacing	2-5-2 mm
Electrode material	Gold (Au)
Shaft diameter	7 F (2.3 mm)
Curve type	see below*
Curve shape	Steerable – FullCircle
Effective catheter length	110 cm
Connection plug	9-pole connector (REDEL), female, 60° coded
Patient cable	PK-109, 111, 112, 142, 143, 147, 150

Patient Cable	Ablation System	Order Number
PK-109	Stockert EP-Shuttle	330079
PK-111	Osypka HAT 300 Smart	330080
PK-112	Medtronic Atakr Medtronic Atakr II	330081
PK-142	Abbott (SJM) IBI-1500 T7 & T11 Abbott Ampere	362442
PK-143	Boston Scientific Maestro 3000	371597
PK-147	BIOTRONIK Qubic RF, Stockert EP-Shuttle	398853
PK-150	Stockert SmartAblate	402668

Curve Type*	Reach R
Red	48 mm
Blue	57 mm
Green	65 mm
Black	73 mm
Cyan	80 mm

Transport and Storage Conditions

Storage	Please store this sterile product in the original sterile packaging
Light	Keep away from sunlight

