# **Extracorporeal Blood Circuit** for Blood Purification Equipment

### Corporate Headquarters

WEIHAI WEIGAO BLOOD PURIFICATION PRODUCTS CO., LTD.

No.20 Xingshan Road, Weihai Torch Hi-tech Science Park, 264210 Weihai, Shandong Province, PEOPLE'S REPUBLIC OF CHINA.

Tel: +86 631 5660598 Fax: +86 631 5660598

#### International Business Headquarters


WEGO HEALTHCARE (SHENZHEN) CO., LTD.

25/F, East Technology Mansion, No.16, Keyuan Road, Nanshan District, Shenzhen, China Tel: +86 755 33892500 Fax:+86 755 33892508 E-mail: info@wego-healthcare.com

www.wego-healthcare.com




## Extracorporeal Blood Circuit for Blood Purification Equipment





### **Product Features**

- Ovarious models with high compatibility: Over ten models are compatible with most dialysis machines in the market.
- A wide range of accessories: Accessories such as infusion line, pre-fill connectors and waste bags are available for you to choose at will, and can meet the different demands.
- O Double lateral branch waste bag: The pre-fill operation is standardized, and the latest SOP requirements and industrial consensus are satisfied.
- The products are certified by CFDA and EU authority.







### Technical specification

Model	Dimension	Features
JRHLS-001	/	Displacement liquid tube for hemofiltration, special pipelines
JRHLL-010	Outer diameter of venous chamber: 20 mm Pumping pipe: 8*12mm	Universal model
JRHLL-020	Outer diameter of venous chamber: 22 mm Pumping pipe: 8*12mm	Universal model
JRHLL-040	Outer diameter of venous chamber: 30 mm Pumping pipe: 8*12mm	For Fresenius large chamber devices