


Stimulating needle Single shot injection

locoplex®

This stimulating needle is insulated over its entire length. Its bevel is completely exposed for a large conduction of current, and for an easier needle insertion.



A choice of 2 exposed bevels to fit your needs
30° bevel needles (short bevel) facilitate perception of the successive subcutaneous layers. It requires a greater force of penetration than 17° bevel needles (long bevel).

Hygiene
A long (60cm) electrical cable to work at distance from sterile field.

Length	17° bevel		30° bevel		Gauge	Packaging
	Code		Code			
25 mm	5194.25		5194.253		23	20 / box
35 mm	5194.35		5194.353		23	20 / box
50 mm	5194.50		5194.503		20/21	20 / box
85 mm	5194.85		5194.853		21	20 / box
100 mm	5194.10		5194.103		21	20 / box
120 mm	5194.12		5194.123		20	20 / box
150 mm	5194.15		5194.153		20	20 / box


VYSET® These references are also available in VYSET, with devices dedicated to ultrasound.
For more details, see page 13.

Catheter with wishbone Continuous injection

techniplex®

Introduced through 17°, 30° or Tuohy needle

techniplex® is introduced through a puncture needle with 17°, 30° or Tuohy bevel. The wishbone allows an excellent control of the catheter's progress through the needle. The technical procedure is therefore rapid and simple with limited manipulations, and is aseptic thanks to the protective sheath.



locoplex®
techniplex® introduction needles have the same advantages as locoplex® needles.

Catheter		18G needle		Needle 17° bevel		Needle 30° bevel		Needle 30° bevel (catheter with stylet and soft closed distal tip (3 lateral eyes))		Needle Tuohy needle (catheter with stylet and soft closed distal tip (3 lateral eyes))		Packaging
Length	G	Length	G									
50 cm	20	35 mm	18	5195.35		5195.353		-		-		10 / box
50 cm	20	50 mm	18	5195.50		5195.503		5195.5031		5195.506		10 / box 5 / box
50 cm	20	100 mm	18	5195.10		5195.103		5195.1031		5195.106		10 / box 5 / box
50 cm	20	120 mm	18	5195.12		5195.123		-		-		10 / box
50 cm	20	150 mm	18	-		-		5195.1531		5195.156		5 / box