

HOLTSCH MED

Quality. Safety. Trust.

## Quickpad

The Quickpad alcohol swab dispenser, being sterile and physiologically verified as harmless, is ideal for cleaning and disinfecting the skin.

The active agent 2 propanol acts as an effective, "mild on skin" disinfecting agent. The "swab tear-off" ready to use system allows for the single use of the swab. This makes Quickpad not only economical but also efficient in its use.

The lid of the Quickpad alcohol swab dispenser container seals air-tight, keeping the alcohol swabs moist and sterile. As a result, the swab dispenser has a particularly long shelf life of 24 months.

Quickpad is supplied ready for use and with its simple operation, can be used by specialists as well as patients for cleaning/disinfecting the skin. An ideal product for diabetics and other users of self-injection syringes. Additionally, it is always possible through the transparent container to see the fill level and it is suitable for self-application by the patient.

Holtsch	T +49 6128 91717-7
Medizinprodukte GmbH	F +49 6128 44742
In den Faltern 13	M info@holtsch-med.com
D-65232 Taunusstein	W holtsch-med.com

