

Knitted vascular prosthesis with collagen Ra K bifurcated CE 1434

KNITTED VASCULAR PROSTHESIS OFFERS:

Direct use

The collagen layer allows the prosthesis to be used without the need for any modification before implantation.

Easy handling

Coloured guiding line for easy orientation in the surgical wound and easy prosthesis handling during the surgery.

Excellent features

The unique design ensures the strength of the prosthesis wall while providing the required porosity, pliability, softness, and ease of sewing.

Simple adjustment

The prosthesis can be cut to the desired shape at the ends with scissors, the edges do not fray and thus allow the creation of a fine anastomosis.

Parameters

Water permeability and leakage $<5 \text{ ml.cm}^{-2}.\text{min}^{-1}$. Impregnation with crosslinked bovine collagen, formaldehyde and glutaraldehyde free.

Indications

The knitted vascular prosthesis with collagen is indicated for surgical treatment of aneurysms and arterial occlusions. The prosthesis is especially suitable for **the reconstruction of the abdominal aorta and peripheral arteries.**

ORDERING INFORMATION (BIFURCATED)

Ø [mm]	Catalogue number length 400 mm*	Catalogue number length 450 mm*	Catalogue number length 500 mm*
6/12	024-0612	024-0612/450	024-0612/500
7/14	024-0714	024-0714/450	024-0714/500
8/16	024-0816	024-0816/450	024-0816/500
9/18	024-0918	024-0918/450	024-0918/500
10/20	024-1020	024-1020/450	024-1020/500
11/22	024-1122	024-1122/450	024-1122/500
12/24	024-1224	024-1224/450	024-1224/500

* These items are not a standard item, produced on the customer's request.

Sterilized by irradiation

Expiration 3 years from the date of manufacture

