



深圳市理邦精密仪器股份有限公司
EDAN INSTRUMENTS,INC.

EDAN\QR\IM70\1.4

PROFORMA INVOICE

To: SRL "Aelo Grup"	Deliver to: SRL "Aelo Grup"
Address: mun. Chisinau MD-2004 str. Petru Movilă 23/5, ap.3	Address: mun. Chisinau MD-2004 str. Petru Movilă 23/5, ap.3
Tel: 37360106007	Tel: 37360106007
Fax: (37322)601491	Fax: (37322)601491
Attn: Rodica Cobzari-Turcan	Attn: Rodica Cobzari-Turcan
E-mail: aelogrup@gmail.com	Invoice No. EDAN-IM-231016FC
From: Shenzhen, China	Date: 16/07/24
To: Moldova	Currency: USD
Payment: T/T	Trade Items: EXW

Part No.	Description of Goods	Quantities	Unit Price(USD)	Amount(USD)
83.60.262543	X12Patient Monitor With Russian Software,CNBP,HDMI, with Sidestream EDAN G2 CO2 function,and standard accessories, with 2IBP function,and standard IBP accessories, with Thermal Recorder module, With User Manual(Russian)	6 pcs		
83.62.560027	SD3 PRO Ultrasonic Pocket Doppler, With 3MHz Waterproof Interchangeable Probe(For FHR Detection), With User Manual(Russian), With Quick Reference(Russian)	10 pcs		
83.61.360647	SE-3(Basic,Colour Screen) Electrocardiograph, With Russian Software, With User Manual(Russian),230V,50Hz	2 pcs		
83.61.360251	SE-601B Electrocardiograph,CE,European standard,EDAN,Russian silicone keyboard, With Russian Software, With User Manual(Russian)	2 pcs		
01.57.471377017	<u>ECG trunk cable,3-</u> <u>lead,Defib,IEC/AHA,2.7m,reusable</u>	6 pcs		
01.57.471069-10	Nellcor SpO2 Extension cable	2 sets		
01.15.30043	Nellcor Reusable Adult SpO2 Sensor(DS- 100A OxiMax) [*CE+FDA+CA+BA*]	4 pcs		
01.57.040437001	Nellcor Reusable SpO2 sensors, Pediatric/Infant 3-40 kg (OXI-P/I)	4 pcs		

01.57.471331019	NIBP Cuff, Larger Adult,34cm-43cm,reusable, With EDAN logo (has to work with EDAN NIBP tube) [*CE+FDA+CA+BA*]	4 pcs
<u>01.57.470330</u>	<u>Adult E9,EDAN,Adult,Reusable,27-35cm, With EDAN Logo (has to work with EDAN NIBP tube) [*CFDA+CE+FDA+CA+BA*]</u>	<u>2 pcs</u>
<u>01.25.471329</u>	<u>NIBP Cuff, Small Adult,20.5cm-28cm,reusable, With EDAN logo (has to work with EDAN NIBP tube)</u>	<u>4 pcs</u>
<u>01.57.471331</u>	<u>NIBP Cuff, Larg Adult,34cm-43cm,reusable, With EDAN logo (has to work with EDAN NIBP tube)</u>	4 pcs
<u>01.57.471328019</u>	<u>NIBP Cuff, Child,16cm-21.5cm,reusable, With EDAN logo (has to work with EDAN NIBP tube)</u>	4 pcs
<u>01.57.471327020</u>	<u>NIBP Cuff, Small Child,13cm-17cm,reusable, With EDAN logo (has to work with EDAN NIBP tube)</u>	4 pcs
<u>01.57.471326018</u>	<u>NIBP Cuff, Infant,10cm-15cm,reusable, With EDAN logo (has to work with EDAN NIBP tube)</u>	
01.57.471971013	Y shape pressure transducer interface cable, BD	2 pcs
01.57.040121011	IBP Pressure transducer kit, BD, disposable (Argon DT-4812 / Merit Medical, 682000)	20 pcs
<u>01.57.07835</u>	<u>EDAN Recording paper</u>	<u>90 pcs</u>
01.57.471275012	Disposable CO2 Sampling line with male luer lock(C0311)	40 pcs
01.57.471286013	Child Nasal CO2 with O2 delivery sampling cannula(C0305)	20 pcs
01.57.471282013	Adult Nasal CO2 sampling cannula(C0307)	20 pcs
01.57.471483017	ECG trunk cable,3-lead,Defib,IEC/AHA,2.7m,reusable	4 pcs
01.57.471069-10	Nellcor SpO2 Extension cable	4 sets
01.15.30043	Nellcor Reusable Adult SpO2 Sensor(DS-100A OxiMax) [*CE+FDA+CA+BA*]	8 pcs
01.57.471331019	NIBP Cuff, Larger Adult,34cm-43cm,reusable, With EDAN logo (has to work with EDAN NIBP tube) [*CE+FDA+CA+BA*]	8 pcs
<u>01.15.0402.53</u>	<u>2-pin Pediatric/Neonate skin temperature probe (2.252K/25°C) SPO2 Finger Sensor,Pediatric,1m,reusable (has to work with EDAN SpO2 Extension Cable for patient monitor)</u>	<u>4 pcs</u>
<u>02.01.210120030</u>	<u>[*CFDA+CE+FDA*]</u>	8 pcs
<u>01.57.471859</u>	<u>ECG Electrodes, child/neonate, disposable, 1 set (50 pcs/set) [*CFDA+CE+FDA+BA*]IBP Pressure transducer kit, BD,</u>	300 pcs
01.57.040121011	disposable (Argon DT-4812 / Merit Medical, 682000)	20 pcs
02.01.210520011	EDAN Water Trap	40 pcs

01.57.471275013	Disposable CO2 Sampling line with male luer lock(C0311)	40 pcs
-----------------	--	--------

01.57.471282013	Adult Nasal CO2 sampling cannula(C0307)	40 pcs
01.57.471377017	ECG cable, 3-lead, clip, Defib, IEC, 3.4m, reusable	2 pcs
01.57.471068020	SpO2 7-pin Extension Cable,2m	1 pc
02.01.210120030	SPO2 Finger Sensor,Adult,1m,reusable (has to work with EDAN SpO2 Extension Cable for patient monitor) [*CFDA+CE+FDA*]	1 pc
02.01.210121025	SpO2 Silicone Soft-tip Sensor,Pediatric,1m,reusable (has to work with EDAN SpO2 Extension Cable for patient monitor) [*CFDA+CE+FDA*]	2 pcs
01.57.471326018	NIBP Cuff, Infant,10cm-15cm,reusable, With EDAN logo (has to work with EDAN NIBP tube)	2 pcs
01.57.471327020	NIBP Cuff, Small Child,13cm-17cm,reusable, With EDAN logo (has to work with EDAN NIBP tube)	2 pcs
01.57.471328019	NIBP Cuff, Child,16cm-21.5cm,reusable, With EDAN logo (has to work with EDAN NIBP tube)	2 pcs
01.57.471329019	NIBP Cuff, Small Adult,20.5cm-28cm,reusable, With EDAN logo (has to work with EDAN NIBP tube)	2 pcs
01.57.471330019	Adult E9,EDAN,Adult,Reusable,27-35cm, With EDAN Logo (has to work with EDAN NIBP tube) [*CFDA+CE+FDA+CA+BA*]	1 pc
01.57.471331019	NIBP Cuff, Larger Adult,34cm-43cm,reusable, With EDAN logo (has to work with EDAN NIBP tube) [*CE+FDA+CA+BA*]	2 pcs
01.54.458134013	Service Manual, English(for X8 X10 X12)	1 pc
01.54.458903012	Service Manual, English(for iX10 iX12 iX15)	2 pcs

Total (EXW) :

Warranty Period	Main Unit 12 months, accessories 12 months.
Payment:	Total prepayment by T/T
Shipment Date:	Within 10-20 working days upon receiving total prepayment.
Remark:	<p>1) Please pay the bank handling charge both in your country bank and Intermediate bank</p> <p>2) Invoice No. & message "Medical Equipment" must be indicated in the Bank form , as per the policy of foreign exchange control of China government</p> <p>3) Please email or fax the bank receipt upon successful payment</p> <p>4) If Buyer choose the EXW trade terms to entrusts a third-party freight forwarding company to pick up the goods at Edan's factory itself or choose the FCA/FOB trade terms to entrusts a third-party freight forwarding company to transfer the cargo rights with EDAN at a designated place, Buyer shall bear all the transportation risk (including the transport accidents, poor seal quality of the transport vehicle, damaged goods and so on) and customs risk of the goods after delivery/transfer, and EDAN shall not bear any responsibility.</p>

Traceability: Buyer shall keep records for each Product resold by Buyer including but not limited to names, address and telephone number of customers, sales records and records of service in order to render it possible to efficiently execute all necessary recalls or mandatory updates of the Products. These records shall be made available to Seller upon Seller's request and maintained by Buyer for eleven (11) years.

Complaint: Buyer shall establish a procedure for handling complaints with respect to the Products as well as maintain a file for each such complaint. Buyer shall inform Seller of all complaints, including the manner in which buyer was able to resolve the issues. Buyer will investigate each complaint with its customer to ascertain such information that will assist Seller in identifying and correcting the complaint. Buyer will provide to Seller written information pertaining to each such complaint wherever possible. Upon request, Buyer shall collect the defective products/parts and send such products/parts to Seller for investigation.

Recall(MDR): If any recall or MDR of any of the Products occurs, Buyer shall assist Seller to recall such Products and deliver such Products to Seller upon Seller's request.

USD A/C Info:

Beneficiary:

EDAN INSTRUMENTS, INC.

Beneficiary Bank:

BANK OF CHINA, SHENZHEN BRANCH

Swift code: BKCHCNBJ45A

USD A/C No: 769257961091

Correspondent Bank:

Bank of China, New York Branch, New York

Swift code: BKCHUS33

ADDRESS: INTERNATIONAL FINANCIAL BUILDING NO.2022,

JIAN SHE ROAD, SHENZHEN, CHINA

Valid Date: 30 days after this pro-forma invoice is issued.



Note: Please remark EDAN **Invoice No.** in remittance Information for wire transfer payment to EDAN account

Authorized Signature

EDAN INSTRUMENTS, INC.

深圳市理邦精密仪器股份有限公司

Tel: 86-755-26882220 Fax: 86-755-26898330 www.edan.com

Add: No.15 Jinhui Rd., Jinsha Community, Kengzi Subdistrict, Pingshan District, Shenzhen, 518122 P.R.China