

Suction catheter

Available in four adult sizes and in two lengths for attachment to endotracheal or tracheostomy tube

Retainable irrigation port cap

Preventing accidental loss of the cap

Soft catheter tip with Murphy eyes

Allows gentle suction of patient secretions

Isolation valve

Isolates catheter for rinsing on the 72-hour systems

When closed, prevents accidental advancement of the catheter

Double swivel elbow

Facilitates the movement of the patient and prevents trauma

Flexible catheter mount

Can be added for increased flexibility and patient comfort

Disconnection wedge

For removal of the closed suction system from the endotracheal and tracheostomy tubes with minimal discomfort to the patient

Catheter length markings

Depth level indicators for easy positioning of the catheter

Improved sleeve

Clear: allowing for constant monitoring of patient secretions
Strong: with improved feel for catheter insertion
Quiet: less noise ensuring the patient is not disturbed

Colour coding

Product size easily identified in line with International Standards

Ergonomic suction valve

Allows perfect control of the applied suction force and a positive locking system to reduce the risk of accidental depression

