

i-gel® supraglottic airway



For use in anaesthesia and resuscitation

i-gel is a truly unique, single use, latex and PVC free second generation supraglottic airway management device that represents the culmination of years of extensive research and development. Made from a medical grade thermoplastic elastomer, i-gel has been designed to create a non-inflatable, anatomical seal of the pharyngeal, laryngeal and perilaryngeal structures whilst avoiding compression trauma.

The i-gel® supraglottic airway is supplied in a fully recyclable protective cradle or cage pack. This unique packaging protects the i-gel in transit and ensures that it maintains its anatomical shape. i-gel is available in seven sizes.

www.i-gel.com

Video available at www.intersurgical.com



15mm connector

Reliable connection to any standard catheter mount or connection


Proximal end of gastric channel

Clearly displayed product information

For quick easy reference. Includes confirmation of size and weight guidance

Position guide (adult sizes only)

Easy confirmation of optimum insertion depth



Gastric channel

The i-gel incorporates a gastric channel (except size 1). It provides an early warning of regurgitation, allows for the passing of a nasogastric tube to empty stomach contents and facilitates venting

Integral bite block

Reduces the possibility of airway channel occlusion

Buccal cavity stabiliser

Aids insertion and eliminates the potential for rotation


Epiglottic rest

Reduces the possibility of epiglottic 'down folding' and airway obstruction

The non-inflatable cuff

Made from a unique soft gel-like material allowing ease of insertion and reduced trauma

Distal end of gastric channel



Code	Description	Size	Weight	Box Qty.
8205000	i-gel, supraglottic airway	5 Large adult	90+kg	25
8204000	i-gel, supraglottic airway	4 Medium adult	50–90kg	25
8203000	i-gel, supraglottic airway	3 Small adult	30–60kg	25
8225000	i-gel, supraglottic airway	2.5 Large paediatric	25–35kg	10
8202000	i-gel, supraglottic airway	2 Small paediatric	10–25kg	10
8215000	i-gel, supraglottic airway	1.5 Infant	5–12kg	10
8201000	i-gel, supraglottic airway	1 Neonate	2–5kg	10



Make an enquiry