

A circular inset image showing a box of microscope slides and several individual slides of different colors (blue, green, yellow, red, orange).

CHONGQING
NEW WORLD TRADING CO., LTD.

PRODUCTS CATALOG

www.cnwtc.com



CONTACT US

ADD: Room10-1, Unit2, No.88, Jianxin East Road,
Wulidian Street, Jiangbei District, Chongqing, China

Tel: 0086 023-67728212

Fax: 0086 023-67729109

E-mail: kathy@cnwtc.com export@cnwtc.com

CHONGQING
NEW WORLD TRADING CO., LTD.

CNWTC promises to offer reliable quality
and professional service forever.



Company Profiles ▼

Chongqing New World Trading Co., Ltd. (CNWTC) was established in 2009, which is a professional manufacturer and exporter in the area of laboratory and medical products. Our factory locates in Haimen City, Jiangsu Province, and only 2 hours' driving to Shanghai port. The geographical advantage allows us to provide convenient and fast delivery.

We supply full range of products which covers laboratory consumables (pipette tip, test tube, centrifuge tube, urine container, cryo tube, sample cup), medical consumables (blood collection tube, swabs, blood lancet, alcohol pad, face mask, disposables gown, nitrile gloves, coverall) and laboratory & medical equipment. And those products were CE and ISO13485 certified.

We have good experience in supplying qualified products and considerate service to customers from all over the world like European, Middle East, Southeast Asia, Africa, South & North America and so on. Looking forward to cooperate with you.



Product Usage


Microcentrifuge tube is a small centrifuge tube, which is used in conjunction with a micro centrifuge for the separation and centrifugation of trace reagents, and is widely used in the field of life sciences and clinical trials.

Features

- Full range volumes from 0.1ml to 10ml, apply to different solutions;
- Different bottom shapes and cap types for various testing purpose;
- Various Colors can be customized for distinguishing different samples;
- Made of high quality polypropylene raw material, enables excellent view of tube contents;
- Clear mold graduation and big writing area makes easy reading and recording (for 0.5ml-10ml);
- Tight and leak-free seal can open easily, even with gloves;
- Withstand high centrifugal forces;
- Bulk pack, RNase DNase DNA Free, EO sterilization are optional.



Various Colors



Model No.	Description	Cap Type	Material	Packing
TYA201	Centrifuge Tube 10ml, Round bottom, mold graduation	Snap Cap	PP	200pcs/bag, 3000pcs/ctn
TYA202	Centrifuge Tube 7ml, Round bottom, mold graduation	Snap Cap	PP	200pcs/bag, 4000pcs/ctn
TYA203	Centrifuge Tube 5ml, Round bottom, mold graduation	Snap Cap	PP	200pcs/bag, 6000pcs/ctn
TYA204	Centrifuge Tube 5ml, Conical bottom, mold graduation	Snap Cap	PP	100pcs/bag, 2000pcs/ctn
TYA205	Centrifuge Tube 2ml, Round bottom, mold graduation	Snap Cap	PP	500pcs/bag, 10000pcs/ctn
TYA206	Centrifuge Tube 2ml, conical bottom mold graduation	Snap Cap	PP	500pcs/bag, 10000pcs/ctn
TYA207	Centrifuge Tube 1.5ml, conical bottom, mold graduation	Snap Cap	PP	500pcs/bag, 10000pcs/ctn
TYA208	Centrifuge Tube 1.5ml, conical bottom, Shield Lid, mold graduation	Snap Cap	PP	500pcs/bag, 10000pcs/ctn
TYA209	Centrifuge Tube 0.5ml, conical bottom, Shield Lid, mold graduation	Snap Cap	PP	1000pcs/bag, 25000pcs/ctn
TYA210	centrifuge tube 0.5ml, conical bottom, mold graduation	Snap Cap	PP	500pcs/bag, 10000pcs/ctn
TYA211	Centrifuge Tube 0.2ml, conical bottom, no Graduation, Flat Lid	Snap Cap	PP	1000pcs/bag, 50000pcs/ctn
TYA212	Centrifuge Tube 0.2ml, conical bottom, no Graduation, convex Lid	Snap Cap	PP	1000pcs/bag, 50000pcs/ctn
TYA213	Centrifuge Tube 0.1ml, conical bottom, no Graduation, Flat Lid	Snap Cap	PP	1000pcs/bag, 50000pcs/ctn



Introduction

Blood collection tubes are used for collection and packaging of blood samples for biochemical and immunological examinations in medical tests, and are used in clinical medical experiments using human venous blood as samples. According to the different objects of blood samples, blood collection tubes are divided into serum blood collection tubes, plasma blood collection tubes and whole blood collection tubes.

Features

- The tube body is made of PET material or glass material, with stable properties;
- The head cap is safe, with good air tightness, small puncture force, and convenient opening and capping;
- The inner wall is extremely smooth and be evenly coated by additive, to ensure the blood samples fully contact after entering the tube;
- Good compatibility, can match the transmission device of various biochemical, immune, blood and other analyzers;
- Labels can be customized according to customer needs.

Category	Test Item	Additive	Item	Cap Color	Specification
Serum blood collection tube	Clinical Biochemistry, Immunology and Serology Test	No additive	Plain Tube	Red	13x75mm: 2ml, 3ml, 4ml, 5ml 13x100mm: 5ml, 6ml, 7ml
		Clot Activator	Clot Activator Tube	Red Orange	16x100mm: 7ml, 8ml, 9ml, 10ml
		Gel&Clot Activator	Gel&Clot Activator Tube	Yellow	13x75mm: 2ml, 3ml, 3.5ml, 4ml 13x100mm: 5ml, 6ml 16x100mm: 7ml, 8ml, 8.5ml
Whole blood collection tube	Haematology Test (Blood Routine Examination)	K2EDTA/ K3EDTA	EDTA Tube	Purple	13x75mm: 2ml, 3ml, 4ml, 5ml 13x100mm: 5ml, 6ml, 7ml 16x100mm: 7ml, 8ml, 9ml, 10ml
	Erythrocyte Sedimentation Rate Test	3.8% Sodium Citrate Buffer (0.129mol/L)	ESR Tube	Black	13x75mm: 1.6ml, 2.4ml, 3.2ml 8x120mm: 1.26ml, 1.6ml, 2.0ml 9x120mm: 1.6ml
Plasma Blood Collection Tube	Coagulation Function Test	3.2% Sodium Citrate Buffer (0.109mol/L)	PT Tube	Blue	13x75mm: 1.8ml, 2.0ml, 2.4ml, 2.7ml, 3ml, 3.6ml, 4ml, 4.5ml 13x100mm: 5ml, 6ml, 7ml 16x100mm: 8ml, 9ml, 10ml Double Wall type: 13x75mm: 1.8ml, 2.7ml
	Clinical Chemistry for Emergency Treatment, Blood Rheology Test	Sodium Heparin/ Lithium Heparin	Heparin Tube	Green	13x75mm: 2ml, 3ml, 4ml, 5ml 13x100mm: 5ml, 6ml, 7ml 16x100mm: 7ml, 8ml, 9ml, 10ml
		Gel&Sodium Heparin/ Gel& Lithium Heparin	Gel&Heparin Tube		13x75mm: 2ml, 3ml, 3.5ml, 4ml 13x100mm: 5ml, 6ml 16x100mm: 7ml, 8ml, 8.5ml
	Glucose and Lactate Test	Sodium Fluoride& Sodium Heparin/ Sodium Fluoride& EDTA/ Sodium Fluoride& Potassium Oxalate	Glucose Tube	Grey	13x75mm: 2ml, 3ml, 4ml, 5ml



Serum Blood Collection Tube

Serum blood collection tubes are designed to provide uncontaminated, undifferentiated raw blood samples for medical testing. After centrifugation, serum can be effectively separated from blood cells and fibrin. Serum tube includes 3 tubes: Plain tube with red cap, Clot Activator tube with red/orange cap, and Gel & Clot Activator tube with yellow cap.



Plain Tube, Red Cap

Plain tubes also called No Additive tube, which is used for collection and packaging of blood samples such as medical testing, biochemistry, immunology, serology, various virus detection, blood bank inspection, etc.

Plain tube can maintain the normal composition of serum for a long time and facilitate the re-examination of serum with good repeatability.



TYK101-203

Model No.	Material	Additive	Size(mm)	Vol	Packing
TYK101	PET	No additive	13x75	2ml, 3ml, 4ml, 5ml	100pcs/rack, 1800pcs/ctn
TYK201	Glass		13x100	5ml, 6ml, 7ml	100pcs/rack, 1800pcs/ctn
TYK102	PET				
TYK202	Glass				
TYK103	PET		16x100	7ml, 8ml, 9ml, 10ml	100pcs/rack, 1200pcs/ctn
TYK203	Glass				



Clot Activator Tube, Red/Orange Cap

Clot Activator tube, also known as pro-coagulant tube, is used to collect serum samples for biochemical and immunological testing, and can be used for rapid outpatient, emergency serum biochemical tests.

The inner wall is extremely smooth and be evenly coated by coagulant.

It has a wide operating temperature range. Blood clots completely within minutes, and centrifugation yields high quality serum.



TYK104-206

Model No.	Material	Additive	Size(mm)	Vol	Packing
TYK104	PET	Clot Activator	13x75	2ml, 3ml, 4ml, 5ml	100pcs/rack, 1800pcs/ctn
TYK204	Glass		13x100	5ml, 6ml, 7ml	100pcs/rack, 1800pcs/ctn
TYK105	PET				
TYK205	Glass				
TYK106	PET		16x100	7ml, 8ml, 9ml, 10ml	100pcs/rack, 1200pcs/ctn
TYK206	Glass				



Gel and Clot Activator Tube, Yellow Cap

Gel & Clot Activator tubes are used for serum biochemistry, drug detection and immunology etc.

Coagulant is finely coated on the inner wall of which can shorten the blood coagulation time. After centrifugation, the serum and blood cells are completely separated to obtain high quality serum.



TYK107-209

Model No.	Material	Additive	Size(mm)	Vol	Packing
TYK107	PET	Gel & Clot Activator	13x75	2ml, 3ml, 3.5ml, 4ml	100pcs/rack, 1800pcs/ctn
TYK207	Glass		13x100	5ml, 6ml	100pcs/rack, 1800pcs/ctn
TYK108	PET				
TYK208	Glass				
TYK109	PET		16x100	7ml, 8ml, 8.5ml	100pcs/rack, 1200pcs/ctn
TYK209	Glass				

Lancing Device

Introduction

The depth of the skin jab is measured by the number of millimeters the needle penetrates the skin. The lower-numbered depth setting is shallower, while higher numbers mean a deeper jab. While a lower-numbered depth setting will be more comfortable for most people, a shallow jab may not work well for people with rough or thick skin. The ideal lancet depth will provide enough blood to saturate a glucose strip without causing long-term skin damage or a great deal of pain.

Features

- Little pain design
- Adjustable 9 penetration depths
- Universal for all standard lancets in the market
- Durable and easy to use

Packing: 1pc/box, 100pcs/ carton

Model No.	TYJ601	TYJ602	TYJ603	TYJ604
Size (Diameter*Length)	15x93mm	15x105mm	16x100mm	15x120mm
Available Colors	Beige/white/black or customize	white/black or customize	Blue/white/black or customize	Black/Beige/white or customize



TYJ601-TYJ604

Blood Collection Needle

Introduction

Blood Collection Needle is a safety-engineered, multi-sample blood collection needle that offers a simple, effective way to collect blood while reducing the possibility of needle stick injuries. Needles are available with different types and different sizes for different blood collection needs.

Pen Type Vacutainer Needles

Material: PP, stainless steel

Sterilization: gamma radiation

Gauge: 18g, 20g, 21g, 22g, 23g

Packing: 100pcs/box, 5000box/carton



TYJ701-TYJ705

Model No.	TYJ701	TYJ702	TYJ703	TYJ704	TYJ705
Size	18G	20G	21G	22G	23G
Diameter	1.2mm	0.9mm	0.8mm	0.7mm	0.6mm
Needle Length	1"or1 1/2" (25mm or 38mm)	1"or1 1/2" (25mm or 38mm)	1"or1 1/2" (25mm or 38mm)	1"or1 1/2" (25mm or 38mm)	1"or1 1/2" (25mm or 38mm)
Color	Pink	Yellow	Green	Black	Blue

Pen Type Flashback Vacutainer Needles

Pen Type Flashback Vacutainer Needle is a multi-sample needle that provides instant confirmation of needle placement, through a unique visual flashback indication.

Material: PP, stainless steel

Sterilization: gamma radiation

Gauge: 21g, 22g, 23g

Packing: 100pcs/box, 5000pcs/carton



TYJ801-TYJ803

Model No.	TYJ801	TYJ802	TYJ803
Size	21G	22G	23G
Diameter	0.8mm	0.7mm	0.6mm
Needle Length	1"or1 1/2"(25mm or 38mm)	1"or1 1/2"(25mm or 38mm)	1"or1 1/2"(25mm or 38mm)
Color	Green	Black	Blue

Butterfly Needles

Material: PP, stainless steel

Sterilization: gamma radiation

Gauge: 20g, 21g, 22g, 23g, 25g

Packing: 100pcs/box, 20box/carton (2000pcs/carton)
or 100pcs/bag, 50bags/carton (5000pcs/carton)



TYJ901-TYJ905

Model No.	TYJ901	TYJ902	TYJ903	TYJ904	TYJ905
Size	20G	21G	22G	23G	25G
Diameter	0.9mm	0.8mm	0.7mm	0.6mm	0.55mm
Needle Length	3/4"or1"	3/4"or1"	3/4"or1"	3/4"or1"	3/4"or1"
Tube Length	18/25/30cm	18/25/30cm	18/25/30cm	18/25/30cm	18/25/30cm
Color	Yellow	Green	Black	Blue	Orange

Butterfly Needle with Holder

Material: PP, stainless steel**Sterilization:** gamma radiation**Gauge:** 20g, 21g, 22g, 23g, 25g**Packing:** 50pcs/box, 1000pcs/carton

TYJ1001-TYJ1005

Model No.	TYJ1001	TYJ1002	TYJ1003	TYJ1004	TYJ1005
Size	20G	21G	22G	23G	25G
Diameter	0.9mm	0.8mm	0.7mm	0.6mm	0.55mm
Needle Length	3/4"or1"	3/4"or1"	3/4"or1"	3/4"or1"	3/4"or1"
Tube Length	18/25/30cm	18/25/30cm	18/25/30cm	18/25/30cm	18/25/30cm
Color	Yellow	Green	Black	Blue	Orange

Safety Type Butterfly Needle



TYJ1101-TYJ1103

Material: PP, stainless steel**Sterilization:** gamma radiation**Gauge:** 21g, 22g, 23G**Packing:** 100pcs/box, 20box/carton (2000pcs/carton) or 100pcs/bag, 50bags/carton (5000pcs/carton)

Model No.	TYJ1101	TYJ1102	TYJ1103
Size	21G	22G	23G
Diameter	0.8mm	0.7mm	0.6mm
Needle Length	3/4"or1"	3/4"or1"	3/4"or1"
Tube Length	30cm	30cm	30cm
Color	Green	Black	Blue



Scalp Vein Set

Introduction

Scalp Vein Set is a specialized device used for drawing blood from a vein and giving medication or intravenous (I.V.) therapy into a vein. It consists of a very fine hypodermic needle, two flexible "wings" -resembling the wings of a butterfly-flexible transparent tubing, and a connector. The connector is attached to another device i.e. to a vacuum tube or collection bag to draw blood, or to extension tubing from an infusion pump or IV bag to deliver medication and fluids.

Normal Type

Material: PP+stainless steel**Sterilization:** gamma radiation**Gauge:** 21g, 22g, 23g**Type:** Luer Lock & Luer Slip**Packing:**

100pcs/box, 20 boxes/carton (2000pcs/carton)

or

100pcs/bag, 50 bags/carton (5000pcs/carton)



TYJ1201-TYJ1206

Model No.	TYJ1201	TYJ1202	TYJ1203	TYJ1204	TYJ1205	TYJ1206
Size	21G	22G	23G	21G	22G	23G
Diameter	0.8mm	0.7mm	0.6mm	0.8mm	0.7mm	0.6mm
Needle Length	3/4"or1"	3/4"or1"	3/4"or1"	3/4"or1"	3/4"or1"	3/4"or1"
Tube Length	18/30cm	18/30cm	18/30cm	18/30cm	18/30cm	18/30cm
Color	Green	Black	Blue	Green	Black	Blue
Connector Type	Luer slip	Luer slip	Luer slip	Luer lock	Luer lock	Luer lock

Safety Type

Material: PP+stainless steel**Sterilization:** gamma radiation**Gauge:** 21g, 22g, 23G**Type:** Luer Lock & Luer Slip**Packing:**

100pcs/box, 20 boxes/carton (2000pcs/carton)

or

100pcs/bag, 50 bags/carton (5000pcs/carton)



TYJ1301-TYJ1306

Model No.	TYJ1301	TYJ1302	TYJ1303	TYJ1304	TYJ1305	TYJ1305
Size	21G	22G	23G	21G	22G	23G
Diameter	0.8mm	0.7mm	0.6mm	0.8mm	0.7mm	0.6mm
Needle Length	3/4"or1"	3/4"or1"	3/4"or1"	3/4"or1"	3/4"or1"	3/4"or1"
Tube Length	18/30cm	18/30cm	18/30cm	18/30cm	18/30cm	18/30cm
Color	Green	Black	Blue	Green	Black	Blue
Connector Type	Luer slip	Luer slip	Luer slip	Luer lock	Luer lock	Luer lock



TYO301



TYO302



TYO303



TYO304



TYO305



TYO306



Product Usage

Sampling swab are used to store microbiological samples and transport them to clinical laboratories.



Features

- Composed by stick, tip and sealing tube;
- The stick material is available in stainless steel, wood, PP, PS or ABS;
- The swab tip material is available in superior absorbent cotton, artificial cotton, polyester fiber;
- The sealing tube material is available in virgin PP/PS;
- Special treatment on the stick makes the tip wrapped on the stick tightly;
- Complete sampling kits for the safe collection, transport and preservation;
- Adapted to sampling from the parts of urethra and so on;
- Used for sampling from the parts of vagina, urogenital canal, etc.

Model No.	Description	Size(mm)	Stick Material	Tip Material	Sterilization	Packing
TYO301	Cotton swab	150	Bamboo/ Wooden/ Plastic	Cotton	Non-sterile	100pcs/bag,20000pcs/ctn
TYO302	Cotton swab	150	Bamboo/ Wooden/ Plastic	Cotton	EO sterile	Individual, 1pcs/bag,10000pcs/ctn
TYO303	Female Swab	12X175	Wooden/Plastic	Cotton	EO sterile	100pcs/bag,2000pcs/ctn
TYO304	Viscose Swab	12x175	PP/PS	Viscose	EO sterile	200pcs/bag,2000pcs/ctn
TYO305	Sponge Swab	125	PP	Sponge	EO sterile	100pcs/bag,10000pcs/ctn
TYO306	Male Swab	12x86	Stainless	Cotton/ Viscose	EO sterile	200pcs/bag, 3000pcs/ctn