







































| Compatible Brand | P/N | Description | Male Connector | Skin/Oral Side |
|-------------------|----------|---|----------------|----------------|
| Air Shield | U6105-1B | Adult Skin surface Temperature probe, 6.86KΩ,3pin,3meter | | |
| | U6105-1D | Pediatric Skin Temperature probe, 6.86KΩ,3pin,3meter | | |
| | U6105-1N | Infant Skin Temperature probe,6.86KΩ,3pin,3meter | | |
| Atom | U6106-1A | Adult Rectal/Esophageal Temperature probe, 2.25KΩ,6pin,3meter | | |
| | U6106-1C | Pediatric Rectal/Esophageal Temperature probe,2.25KΩ,6pin,3meter | | |
| | U6106-1M | Infant Rectal/Esophageal Temperature probe,2.25KΩ,6pin,3meter | | |
| | U6106-1B | Adult Skin surface Temperature probe,2.25KΩ,6pin,3meter | | |
| | U6106-1D | Pediatric Skin Temperature probe,2.25KΩ,6pin,3meter | | |
| | U6106-1N | Infant Skin Temperature probe,2.25KΩ,6pin,3meter | | |
| Artema | U632-1A | Adult Rectal/Esophageal Temperature probe, 2.25KΩ,10pin,3meter | | |
| | U632-1C | Pediatric Rectal/Esophageal Temperature probe,2.25KΩ,10pin,3meter | | |
| | U632-1M | Infant Rectal/Esophageal Temperature probe,2.25KΩ,10pin,3meter | | |
| | U632-1B | Adult Skin surface Temperature probe,2.25KΩ,10pin,3meter | | |

A~A

| Compatible Brand | P/N | Description | Male Connector | Skin/Oral Side |
|------------------|---------|--|----------------|----------------|
| Artema | U632-1D | Pediatric Skin Temperature probe,2.25KΩ,10pin,3meter | | |
| | U632-1N | Infant Skin Temperature probe,2.25KΩ,10pin,3meter | | |
| Biolight | U609-1A | Adult Rectal/Esophageal Temperature probe, 2.25KΩ,2pin,3meter | | |
| | U609-1C | Pediatric Rectal/Esophageal Temperature probe,2.25KΩ,2pin,3meter | | |
| | U609-1M | Infant Rectal/Esophageal Temperature probe,2.25KΩ,2pin,3meter | | |
| | U609-1B | Adult Skin surface Temperature probe,2.25KΩ,2pin,3meter | | |
| | U609-1D | Pediatric Skin Temperature probe,2.25KΩ,2pin,3meter | | |
| | U609-1N | Infant Skin Temperature probe,2.25KΩ,2pin,3meter | | |
| | U609-1F | Biolight Temperature interface cable, Rectangle 2Pin,3meter | | |
| Comen | U638-1A | Adult Rectal/Esophageal Temperature probe,10KΩ,2pin,3meter | | |
| | U638-1C | Pediatric Rectal/Esophageal Temperature probe,10KΩ,2pin,3meter | | |
| | U638-1M | Infant Rectal/Esophageal Temperature probe, 10KΩ,2pin,3meter | | |
| | U638-1B | Adult Skin surface Temperature probe,10KΩ,2pin,3meter | | |

A~C

| Compatible Brand | P/N | Description | Male Connector | Skin/Oral Side |
|------------------|---------|--|---|---|
| Comen | U638-1D | Pediatric Skin Temperature probe,10KΩ,2pin,3meter |  |  |
| | U638-1N | Infant Skin Temperature probe,10KΩ,2pin,3meter |  |  |
| Creative | U610-1A | Adult Rectal/Esophageal Temperature probe,10KΩ,2pin,3meter |  |  |
| | U610-1C | Pediatric Rectal/Esophageal Temperature probe,10KΩ,2pin,3meter | |  |
| | U610-1M | Infant Rectal/Esophageal Temperature probe,10KΩ,2pin,3meter | |  |
| | U610-1B | Adult Skin surface Temperature probe,10KΩ,2pin,3meter | |  |
| | U610-1D | Pediatric Skin Temperature probe,10KΩ,2pin,3meter | |  |
| | U610-1N | Infant Skin Temperature probe,10KΩ,2pin,3meter | |  |
| Datascope | U614-1A | Datascope Adult Rectal/Esophageal Temperature probe,2.25KΩ,5pin,3meter |  |  |
| | U614-1C | Rectal/Esophageal Temperature probe,2.25KΩ,5pin,3meter | |  |
| | U614-1M | Infant Rectal/Esophageal Temperature probe,2.25KΩ,5pin,3meter | |  |
| | U614-1B | Adult Skin surface Temperature probe,2.25KΩ,5pin,3meter | |  |
| | U614-1D | Pediatric Skin Temperature probe,2.25KΩ,5pin,3meter | |  |

| Compatible Brand | P/N | Description | Male Connector | Skin/Oral Side | |
|--------------------------|--|---|---|---|---|
| Datascope | U614-1N | Infant Skin Temperature probe,2.25KΩ,5pin,3meter |  |  | |
| Drager/Air Shields C2000 | U615-2N | Drager Infant Skin Temperature Probe,Mini5Pin,9.1KΩ,3meter |  |  | |
| GE | U602-1A | Adult Rectal/Esophageal Temperature probe,2.25KΩ,11pin,3meter |  |  | |
| | U602-1C | Pediatric Rectal/Esophageal Temperature probe,2.25KΩ,11pin,3meter | |  | |
| | U602-1M | Infant Rectal/Esophageal Temperature probe,2.25KΩ,11pin,3meter | |  | |
| | U602-1B | Adult Skin surface Temperature probe,2.25KΩ,11pin,3meter | |  | |
| | U602-1D | Pediatric Skin Temperature probe,2.25KΩ,11pin,3meter | |  | |
| | U602-1N | Infant Skin Temperature probe,2.25KΩ,11pin,3meter | |  | |
| | U602-2C | GE V100 Oral Temperature Probe,Crystal head Connector | |  |  |
| | U602-1F | GE Temperature interface cable,Rectangle 2Pin,3meter | |  |  |
| | U602-1G | GE Temperature Interface Cable,3.5mm Female | |  |  |
| | U602-1J | GE Temperature extension cable,to YSI 400 series & YSI 700 series | |  |  |
| U602-3J | GE Temperature interface cable,6.3mm Female,0.5meter |  |  | | |