Arteriofix

Arterial catheter kits

	Arteriofi	x					
	Arterial catheter kit for invasive monitoring of blood pressure and blood sampling using the Seldinger technique						
	 Introducer needle made of stainless steel smoothly tapered hub transition for easy insertion of guide wire integrated design, providing smooth transition from guide wire to needle and needle to catheter to ensure minimal trauma at puncture site Luer lock fitting 						
	 Catheter made of FEP for pressure transmission and good physical and mechanical properties of the material excellent tissue and blood compatibility smooth surface prevents thrombogenesis special tip for easy insertion FEP catheters have good gliding characteristics during insertion fixation wings at proximal end colour-coded Luer lock-fittings 						
	 Guide wire made of stainless steel – flexible straight tip – designed according to the different inner diameters of catheters 						
	colour	size	capillary: mm [Ø]; length [mm]	needle: Ø [G; mm]; length [mm]	wire: length [cm] – Ø [Inch]	Code- Number	Sales unit-pcs.
	Ateriofix Catheter kits						
	yellow	22G/ 80 mm	0.87; 80	22; 0.80 x 50	25-0.021	<mark>5206316</mark>	20
	red	20G/ 80 mm	1.08; 80	20; 0.95 x 50	25-0.025	5206324	20
	white	20G/160 mm	1.08; 160	20; 0.95 x 70	40-0.025	5206332	20

1.29; 160

18; 1.30 x 70

40-0.035

5206359

20

18G/160 mm

grey