

| DATE | 28.03.2025 | | | | |
|----------|-----------------|--|--|--|--|
| DOC.NO | FRM-KLT-BYM-069 | | | | |
| PAGE NO | 1 | | | | |
| REV.NO | | | | | |
| REV.DATE | | | | | |

| TECHNICAL DATA SHEET | | | | | | | | | |
|--|-------------------------|---|--|---|---------------------|-------------------------|------------|------------|--|
| PRODUCT: Surgical set for laborious interventions (long duration), sterile | | | | | | | | | |
| | | | | | | | | | |
| | | | 72 Surgical set for laborious interventions (long duration), sterile | | | | | | |
| | | | 1 Instrument table field with SMS reinforced area 150x190cm (+/-10cm) - 1 pc. | | | | 150x190 cm | | |
| Description of Product: | | 2 Mayo table field with SMS reinforcement 80x140cm (+/-10cm) - 1 pc. | | | | 80x140 cm | - | | |
| | | 3 Absorbent napkin - not to leave lint - 30x35cm (+/- 10cm) - 2 pcs. U-shaped field. 200x310cm (+/-10cm) - 1 pc. • Incision area - 10x110cm (+/-5cm); • positioning tape - on | | | | 30x35 cm | + | | |
| | | 4 the entire surface of the incision area, impregnated in fabric, width >5cm; • Absorbent area - on the entire surface of the incision, resistant to liquid penetration of min. 150cm of H20/cm2 | | | | 200x310 cm | | | |
| | | 5 positioning tape 10*50 (+/-5cm) - 1 pc | | | | 10x50 cm |] | | |
| | | Lateral field. 95*105cm (+/- 10cm) - 2 pcs • positioning tape - on the entire surface of the incision area, impregnated in fabric, width > 5cm; • Absorbent area - on the entire surface of the incision; resistant to liquid penetration of min 150cm H20/cm2 • SMS material on the entire surface of the field, with liquid resistance min. 150cm H20/cm2 | | | | 95x105 cm | | | |
| | Raw Mater | ials: | SMS /SMS + Reinforcement / Cellulose | | | | | | |
| | Product Col | our: | Medical Blue | | | | | | |
| | Reference C | ode: | SP-03001-67 | | | | | | |
| Weight in Grams: | | | 40 gsm SMS / 35 gsm SMS + Reinforcement / 60 gsm (Cellulose) | | | | | | |
| | Package | : | Flat Pouch | , | | | | | |
| | Pr | roduct: Medical Devi | ce Directive MDR 2017/745 EC | CC ANNEX XI A ALL VERS Product Shelf Life - 5 Year | | ials used do not contai | in metal. | | |
| Product Materials | | | | | PROPERTIES | | | | |
| Unit / Size | | 1 | 2 | 3 | 4 | | | | |
| 1 | SMS | 150x190 cm | _ | | | | | | |
| 1 | SMS | 80x140 cm | | | | | | | |
| 2 | Cellulose | 30x35 cm | 5 | 6 | 7 | | | | |
| 1 | SMS | 200x310 cm | | - | / | | | | |
| 1 | *** | 10x50 cm | | | | | | | |
| 2 | SMS + Reinforecement | 95x105 cm | | | | | | | |
| | lerances: | +/- 2% cm | | | Package Information | | | | |
| Me | asurement: | cm | | | | | | | |
| | | | The products in the sterilization bag are double-packaged to reduce all risks during transportation. Double packaged products are put into | | | | | e put into | |
| | | | Baymed's standard sized carton; dimensions are as follows: | | | | | | |
| | | | Height = 44 cm ; Length = 40 cm ve Width = 60 cm . | | | | | | |
| | | | | | | | | | |
| Preparation Date | | | QUALITY CONTROL APPROVAL | | | | | | |
| | | | | | | | | | |