



**Suction Catheter** | Code : 2001  
*(Plain Connector)*

• **Description**

• Suction Catheter is used to provide access to respiratory tract to remove secretion with maximum comfort to the patient.

• **Features**

- 1- For removal of Secretion from trachea and bronchial region.
- 2- Soft, frosted and kink resistant medical PVC tubing.
- 3- Atraumatic, soft & rounded open tip with two lateral eyes.
- 4- Color coded connector for easy identification of size.
- 5- Sizes: 6 FG to 20 FG.

Code	Size	Color
2001-6	6 FG	Light Green
2001-8	8 FG	Blue
2001-10	10 FG	Black
2001-12	12 FG	White
2001-14	14 FG	Green
2001-16	16 FG	Orange
2001-18	18 FG	Red
2001-20	20 FG	Yellow

**Suction Catheter** | Code : 2002  
*(Thumb Control Connector)*

• **Description**

• Suction Catheter is used to provide access to respiratory tract to remove secretion with maximum comfort to the patient the connector regulates the vacuum by covering or partially covering a vent port with the thumb of the clinician who is doing the suctioning of the patient.

• **Features**

- 1- For removal of Secretion from trachea and bronchial region.
- 2- Soft, frosted and kink resistant medical PVC tubing.
- 3- Atraumatic, soft & rounded open tip with two lateral eyes.
- 4- Color coded connector for easy identification of size.
- 5- Sizes: 5 FG to 20 FG.

Code	Size	Color
2002-5	5 FG	Grey
2002-6	6 FG	Light Green
2002-8	8 FG	Blue
2002-10	10 FG	Black
2002-12	12 FG	White
2002-14	14 FG	Green
2002-16	16 FG	Orange
2002-18	18 FG	Red
2002-20	20 FG	Yellow