

## FOSUN TRADE

Provide Global Customers with  
One-stop Integrated Trade Solution



Bund Finance Center, 600 Zhongshan No.2 Road(E), Shanghai 200010, P.R.China

Tel +86 21 23156666

Fax +86 21 61303288

[www.fosuntrade.com](http://www.fosuntrade.com)

## Comprehensive Health



# CONTENTS

## 01 | About Fosun Trade

## 02 | PPE

- Gloves
- Medical Face Mask
- Protective Clothing
- Isolation / CPE Gown
- Face Shield / Caps
- Scrub Suit / Shoe Cover

## 03 | In Vitro Diagnostics

- Qualitative Rapid Test
- Cancer Screening
- Molecular Diagnosis
- POCT
- Fertility (LH & HCG & FSH)

### Veterinary Product

- IVD Diagnosis
- Life Support Equipment
- Dental Equipment
- Infusion Pump
- Imaging Equipment

## 04 | Life Sciences & Lab equipment and Consumable

- Lab Equipments
- Food Quick Testing
- Molecular biology consumables
- Cell Biology Consumables
- Microbiology Consumables
- Pathology Consumables
- Vaccine / Syringe & Needle
- Oral & Enteral / Diabetes Care
- Blood Collection / Infusion Therapy
- First Aid Kit
- Medical Dressing

## 05 | Family Consumer

- Home Care
- Oxygen Concentrator/  
UV-Ozone Disinfection Cart

# SPECIALIZED GLOBAL TRADE PLATFORM OF FOSUN

Provide global customers with One-stop Integrated Trade Solution

## About FOSUN TRADE

Fosun Trade is committed to Fosun's mission of "creating happier lives for families worldwide", and integrates global high-quality supplier resources and provides one-stop integrated trade solutions for global customers to create value on demand, becoming the builder of global family happiness business

## MISSION

Integrate global high-quality supplier resources and provide one-stop integrated trade solutions for global customers to create value on demand

## VISION

Rooted in China, serving one billion family customers worldwide, and creating a healthy, happy and affluent happiness ecosystem

Who  
We Are

Global Customers **4000+**

Supply PPE Material **1301+**

Country Area Covered **30+**

# ABOUT FOSUN TRADE

# Personal Protective Equipment



## PPE

Personal Protective Equipment

### ■ Gloves

#### Disposable Nitrile Exam Gloves

CE Foreal

-Material: Nitrile Rubber (NBR)  
-Feature: Ambidextrous, latex-free, beaded cuff, excellent softness and wearing fitness  
-Standards: FDA, CE, EN 455, food contact, PPE CAT III



FTPPEG-001

size	color
XS,S,M,L,XL	Blue, black, white, violet, etc.

#### Disposable Vinyl Exam Gloves

CE Foreal

-Material: PVC  
-Feature: Ambidextrous, latex-free, beaded cuff, DOP or DEHP free, excellent softness and wearing fitness  
-Standards: FDA, CE, EN 455, food contact, PPE CAT III



FTPPEG-002

size	color
XS,S,M,L,XL	Clear, blue, black, etc.

#### Disposable Latex Exam Gloves

CE Foreal

-Material: Natural Rubber  
-Feature: Ambidextrous, beaded cuff low extractable proteins residue  
-Standards: FDA, CE, EN 455, food contact, PPE CAT III



FTPPEG-003

size	color
XS,S,M,L,XL	Natural color

#### Latex Sterile Surgical Gloves

CE Foreal

-Material: Natural Rubber  
-Feature: Sterilization, anti-slip, packaged in pairs, dedicated left/ right  
-Standards: FDA, CE



FTPPEG-004

size	color
6, 6.5, 7, 7.5, 8, 8.5, 9	Natural color

### ■ Medical Face Mask

#### Medical Protective Masks

Medical Masks - Surgical Masks - Surgical N95 Respirators - Face Mask Earloop Type - Face Mask Earloop Type II/IR - Level 1/2/3 Medical Surgical Mask



FTPPEMFM-001

Classification
general medical masks, surgical masks medical protective masks

### ■ Protective Clothing

#### Disposable Protective Clothing/Lab Coat

Medical purpose/examination, industrial purpose/ PPE laboratory, healthcare and nursing, housekeeping beauty industry, IT industry, food processing



FTPPEPC-001

size	color
S,M,L,XL,XXL	White

### ■ Isolation/CPE Gown

#### Isolation Gown

CE Foreal

Medical purpose/ Examination, Industrial purpose/ PPE, Laboratory, Healthcare and nursing



FTPPEIG-001

size	color
S,M,L,XL,XXL	Blue, yellow, green etc

#### CPE Gown

CE Foreal

Medical purpose/examination, industrial purpose/PPE laboratory, healthcare and nursing, beauty industry, IT industry, food processing



FTPPECPEG-001

size	color
S,M,L,XL,XXL	Blue

## ■ Face Shield

### Face Shield

Anti - scratch  
anti - impact  
anti - dust  
anti - splash



FTPPEFS-001

Packing specifications	Carton size
10PCS/BAG,20BAG/CTN	65x35x30cm

## ■ Caps

### Caps

CE Foreal

Boffant caps are made of multi-ply polypropylene material double-stitched elastic band keeps caps flat and comfortable on the forehead available in a variety of colors and sizes



FTPPEC-001

Material	Type
10g/20g SPP	Boffant caps/Hoods/Surgical Caps

## ■ Scrub Suit

### Scrub suit

Bade with comfortable and breathable nonwoven that's made to last well through long hours of being up in the OT or down in the ward. reliable protection, safe, and fit for use



FTPPESS-001

Material	Type
35g/43g SMS	Green,Blue,Purple

## ■ Shoe Cover

### Shoe Cover

CE Foreal

Provide high-level fluid protection for use in vivariums, pharmaceutical research and manufacturing medium-weight coated polypropylene material attached ties at ankle for secure fit seamless bottom



FTPPESSC-001

Material	Type
Breathable film/SPP/PE	Long/short

In  
Vitro  
Diagnostics



Qualitative Rapid Test

HIV Antibody/Antigen Rapid Test

FTIVDQRT-001

CE (self-testing)

For qualitative detection of HIV-1, HIV-2, O subtype HIVp24 antigens or antibodies (IgG, IgM and IgA) in whole blood serum or plasma, to assist in the diagnosis of HIV infection



Package	Time to Results
25T/Kit, 50T/Kit	10~15minutes

Influenza Virus (Flu) Rapid Test

FTIVDQRT-002

FDA CE (self-testing)

Test item: FluA, FluB  
Sample type: throat, nasal swab nasopharyngeal aspirate  
Storage: room temperature, 24 months



Package	Time to Result
20T/Kit	10~15 minutes

Dengue Rapid Test

FTIVDQRT-003

CE (self-testing)

Test item: Dengue NS1, Dengue IgM, Dengue IgG  
Sample type: whole blood, serum, plasr  
Storage: room temperature, 24 months



Package	Time to Result
25T/Kit, 40T/Kit	15~20 minutes

Helicobacter Pylori (HP) Self-test

FTIVDQRT-004

CE (self-testing)

Accurate: as accurate as the test results of professional medical institutions  
Fast: only 3 steps and results in 10 minutes  
Simple: no need to take medicine before the test



Features	Time to Result
Painless, non-invasive, home-based	10 minutes

Fecal Occult Blood (FOB) Self-test

FTIVDQRT-005

CE (self-testing)

Bad habits (often eating smoked or soured meat)  
Diabetic patients  
Those with a history of enteropathy  
Family history of bowel cancer (76% risk for direct relat)



Features	Time to Result
Painless, non-invasive, home-based	10 minutes

Cancer Screening

Alpha-Fetoprotein (AFP) Rapid Test

FTIVDCS-001

CE (professional)

The Alpha-Fetoprotein Rapid Test is a rapid chromatographic immunoassay for the qualitative detection of AFP in serum, plasma or whole blood to aid in the diagnosis of Hepatocellular Carcinoma or fetal open neural tube defects



Package	Format
25T/Kit, 50T/Kit	strip, cassette

ProCalcitonin Antigen (PCT) Rapid Test

FTIVDCS-002

CE (professional)

The Procalcitonin Rapid Test is used for semi-quantitative determination and monitoring of PCT concentrations in serum, plasma or whole blood specimens



Package	Format
25T/Kit, 50T/Kit	strip, cassette

Carcinoembryonic Antigen (CEA) Rapid Test

FTIVDCS-003

CE (professional)

The Carcinoembryonic Antigen Rapid Test is intended as an aid in the assessment and diagnosis of colorectal cancer and other types of cancers, such as pancreatic, stomach breast lung, and certain types of thyroid and ovarian cancer



Package	Format
25T/Kit, 50T/Kit	strip, cassette

Prostate Specific Antigen (PSA) Rapid Test

FTIVDCS-004

CE (professional) (self-testing)

The Prostate Specific Antigen Rapid Test is intended for use at point of care facilities to measure total prostate specific antigen (tPSA)



Package	Format
5T/Kit, 10T/Kit, 20T/Kit	strip, cassette

Molecular Diagnosis

Automatic Nucleic Acid Extractor

FTIVDMD-001

CE

Display screen: 7/10 inch colour Android screen  
Operation interface: Android system  
Magnetic bead recovery rate: 98%  
Purification space difference CV<3%  
Data storage >500 groups of programs



Purification sensitivity
The positive detection rate of 100 copies of samples>95%

Real-Time Fluorescence Quantitative PCR

FTIVDMD-002

CE

Dimensions: 525mmx320mmx420mm(L\*W\*H)  
Capacity: 96x0.2ml tubes, 12x8-tube strips  
1x96-well plate  
Format: 96-well plate, 8-tube strips, 0.2 ml tubes  
Reaction Volume Range: 10~100µL



Weight	Accuracy
27KG	≤3%

Safety Goggles

POCT Immunofluorescence Analyzer (Portable)

FTIVDSG-001

Samples: serum, plasma, whole blood nasal swab  
Screen: 7 inch touch screen  
Data storage: 10000  
Battery:lithium battery  
Communication: RS232 (LIS/HIS), USB



size	Weight
225mmx152mmx105mm	1.3Kg

POCT Immunofluorescence Analyzer (8 Channels)

FTIVDSG-002

Samples: serum, plasma, whole blood nasal swab  
Screen: 10.1 Inch hd touch screen  
Data storage: 10000  
Test speed: 120t/h  
Communication: RS232 (LIS/HIS), USBWifi, Bluetooth



Dimensions	Weight
390mmx290mmx 205mm	7.5Kg

POCT Molecular Tester

FTIVDSG-003

CE

Easy to use: 1 minute for sample and cartridge loading  
Multiplexed: 1 test with up to 50 targets\*  
Safe: Fully enclosed cartridges minimize the risk of contamination



Fast
Sample reports available within an hour

POCT Immunofluorescence Analyzer (Handhold)

FTIVDSG-004

CE

Samples: serum, plasma, whole bloodnasal swab  
Screen: 5.0 Inch touch screen  
Data storage: 50000  
Battery:lithium battery (super standby time)  
Communication: Type-C, 4G, USB, Wifi



Size	Weight
191 mmx84.5mmx45mm	600g

**Fertility(LH&HCG&FSH)**

**POCT Real-Time PCR**

FTIVDSG-005



·Dimensions: 300mm\*260mm\*180mm(L\*W\*H)  
·Screen: LED 7 inch color touch screen  
·Capacity: 8/16x0.2ml (5-100µl), 8 strips tubes  
·Test dynamic range: 10-10<sup>10</sup> Copies

Uniformity	Heating rate
≤0.1 C	≥6.5 C/s (Max)

**Luteinizing Hormone (LH) Rapid Test**

FTIVDF-001



It is used to detect luteinizing hormone (LH) levels in women's urine to predict the time of ovulation, to guide women to choose the timing of conception or to guide contraception during the safe period

Format	Package
strip, cassette, midstream	1T/Kit, 2T/Kit, 5T/Kit, 10T/Kit, 50T/Kit

**Human Chorionic Gonadotropin (HCG) Rapid Test**

FTIVDF-002



It is used to detect the content and changes of HCG in the urine of women of childbearing age and to determine the auxiliary diagnosis of early pregnancy

Format	Package
strip, cassette, midstream	1T/Kit, 2T/Kit, 5T/Kit, 10T/Kit, 50T/Kit

**Follicle-Stimulating Hormone (FSH) Rapid Test**

FTIVDF-003



To provide useful information to women preparing for pregnancy or concerned about ovarian health, function and fertility, and to aid in the diagnosis and treatment of pituitary and gonadal disorders

Format	Package
strip, cassette, midstream	1T/Kit, 2T/Kit, 5T/Kit, 10T/Kit, 50T/Kit

**Veterinary Auto Hematology Analyzer**

FTIVDVP-005



10.4 inches TFT touch screen, Support up to 20 customized animal detection modes, Build-in thermal printer and USB external printers, Small package reagents

**Cow Blood Ketone Analyzer**

FTIVDVP-006



This device is compact and used to detect the blood ketone content of cows, with results in 9 seconds

**Life Support Equipment**

**Veterinary Monitor**

FTIVDVP-007



Multi-Parameter Veterinary Monitor, High resolution color TFT-LCD Screen, Defibrillation and electrosurgical interference protection

**Veterinary Mini Monitor**

FTIVDVP-008



The device is compact, covering oxygen saturation, pulse rate, temperature, blood pressure and other parameters, and is suitable for monitoring the vital signs of small animals such as dogs or cats

**Veterinary Product**

**Veterinary Immunofluorescence Quantitative Analyzer**

FTIVDVP-001



Using fluorescence immunochromatography technology, this device quantitatively detects common diagnostic indicators such as inflammation, endocrine, reproductive, kidney disease, and heart disease in animals, and produces results within 3-15 minutes

**Veterinary Rapid Test**

FTIVDVP-002



Based on colloidal gold immunochromatography technology, this kit assists in the diagnosis of infectious diseases, parasites, and other infectious diseases in animals, and produces results within 10 minutes

**Veterinary Microfluidic Coagulation Analyzer (electrochemistry)**

FTIVDVP-003



The use of microfluidic electrochemical methodology to detect coagulation factors (APTT, PT, TT, FIB, ACT) in animals has reference significance for the evaluation of coagulation function in dogs and cats

**Veterinary Chemistry Analyzer**

FTIVDVP-004



The equipment adopts Microfluidics technology, and the reagent items include kidney function, liver function, electrolyte, health examination set

**Veterinary ICU Cage**

FTIVDVP-009



The combination of traditional hospital inpatient cages and Multifunctional oxygen treatment system. Application: Emergency Care, Cardiopulmonary Disease, Newborn Care, Postoperative Recovery, Isolation of Infectious Diseases, Elderly Animal Care

**Veterinary Anesthesia Ventilator**

FTIVDVP-010



It is a high-end anesthesia machine used in animal hospitals. The design of the product adheres to the concept of anesthesiologists-oriented, completely abandoning the design ideas of traditional animal anesthesia machines that only considers product costs

**Veterinary Respirator**

FTIVDVP-011



Animal anesthesia ventilator is suitable for automatic respiratory control, respiratory management and respiratory monitoring of different animals during surgery

**Veterinary Blood Pressure Monitor**

FTIVDVP-012



The Blood Pressure Monitor is intended to be used on companion animals, specifically dogs and cats, to obtain oscillometric blood pressure measurements, including systolic, diastolic, mean arterial pressure, and heart beats per minute

## Dental Equipment

### Veterinary Mobile Dental X-ray

Sunnisi

Built-in imported bulb tube, Green X-ray machine, no protection required, High frequency DC voltage, Ergonomic design, Self-contained memory function, Compatible with strong adaptation, Color LCD display



FTIVDVP-013

### Veterinary Dental Scanner

Sunnisi

Support X-ray Photography and scanning. Provides high-resolution images for diagnosis and faster and more accurate diagnosis



FTIVDVP-014

### Veterinary Treatment Unit

Sunnisi

Stable power, Stable performance, Excellent silence, Knob, Hanger



FTIVDVP-015

## Infusion Pump

### Veterinary Syringe Pump

Sunnisi

The veterinary syringe pump is double CPU design. It is open system and can accept different syringe brand after calibration. The software compensation error is within 2% after calibration



FTIVDVP-016

### Veterinary Infusion Pump

Sunnisi

Veterinary infusion pump is with the classic vertical design. The software compensation error is within 5% after calibration



FTIVDVP-017

### Veterinary Blood Infusion Warmer

Sunnisi

Veterinary Blood and Infusion Warmer is to warm the blood or fluid and to relieve the vets hyperthermia; it can transmit heat to IV set to reach the preset temperature through the heating board; deviation is within 1°C; it is widely used for the vet blood or infusion transfusion



FTIVDVP-018

### Veterinary Warming System

Sunnisi

Veterinary Warming System utilizes the convection of air to transmit the hot air generated by the main unit to the through the pipeline, so that the heat energy is evenly distributed and transmitted to the blanket to achieve heat preservation or heating effect



FTIVDVP-019

## Imaging Equipment

### Veterinary Endoscope

Sunnisi

HD image, the details of mucosa can be observed more clearly. Supply workstation software can print report directly. Snap and record, RAM 64GB. φ 8mm insertion tube, φ 2.8mm Working channel



FTIVDVP-020

### Veterinary X-ray System

Sunnisi

Product compose with High Frequency X-ray system, software, flat panel



FTIVDVP-021

### Veterinary Ultrasound

Sunnisi

Handheld ultrasound, compact and portable, matched with different probes, can be used for ultrasound diagnosis of dogs, cats, cows, sheep, and horses. Leading digital beam former technology, CNC Dynamic Frequency Scanning Technology



FTIVDVP-022




# Life Sciences & Lab equipment and Consumable

## Life Sciences & Lab Equipment and Consumable

### ■ Lab Equipments

#### Atomic Absorption Spectroscopy (AAS)

- Atomic Absorption Spectroscopy (AAS)
- Inductively Coupled Plasma Emission Spectrometer (ICP-OES)
- Inductively Coupled Plasma Mass Spectrometer (ICP-MS)



FTLSLEACLB-001 | FTLSLEACLB-002 | FTLSLEACLB-003



#### Chromatography Analysis

- Gas Chromatograph
- liquid chromatography



FTLSLEACLB-004 | FTLSLEACLB-005



#### Mass Spectrometry Analysis

- Gas Chromatography-Mass Spectrometry (GC-MS)
- Liquid Chromatography-Mass Spectrometry (LC-MS)



FTLSLEACLB-006 | FTLSLEACLB-007



#### Molecular Spectroscopy Analysis

- Ultraviolet-Visible Spectrophotometer (UV-Vis)
- Infrared Spectrometer



FTLSLEACLB-008 | FTLSLEACLB-009



## Thermal Analysis and Elemental Analysis

- Differential Scanning Calorimeter (DSC)
- Thermogravimetric Analyzer (TGA)
- Elemental Analyzer



FTLSLEACLB-010 | FTLSLEACLB-011 | FTLSLEACLB-012

## Sample Pretreatment

- Microwave Digestion System
- Solid-Phase Extraction



FTLSLEACLB-013 | FTLSLEACLB-014

## Sample Storage and Handling

- Medical Refrigerator and Freezer
- Biological Safety Cabinet
- CO2 Incubator and Constant Temperature Incubator



FTLSLEACLB-015 | FTLSLEACLB-016 | FTLSLEACLB-017

## Microscope

- Compound High-Power Microscope
- Stereo Microscope
- Digital Microscope



FTLSLEACLB-018 | FTLSLEACLB-019 | FTLSLEACLB-020

## Food Quick Testing

### Colloidal Gold Rapid Test Card

Can be use for detecting biological toxin, heavy metal, drug residue, forbidden addition and easily-abused food additive



FTLSLEACFQT-001

Model	Unit
48T, 96T	box

### Microbial Test Strips

Can be use for detecting microorganism



FTLSLEACFQT-002

Model	Unit
25T	package

### ELISA Test Kit

Can be use for detecting microorganism



FTLSLEACFQT-003

Model	Unit
25T, 50T	box

### Physicochemical Test Reagents

Can be use for detecting biological toxin and nutritional component



FTLSLEACFQT-004


Model	Unit
48T, 96T	box

## Molecular Biology Consumables

### Pipette Tips

FDA CE Foptima

Pipette tips are biological consumables specially designed for quantitative pipetting operations with pipettes. According to the classification of pipettes and different experimental scenes, pipette tips are divided into standard tips of different ranges, extended tips, and low retention tips.




FTLSLEACMBC-001

Specification	Character
10uL, 20uL, 50uL, 100uL, 200uL, 300uL, 1000uL, 1250uL, 5mL, 10mL	Non-filter, Filter, bulk pack, rack pack

### PCR Tubes/Plate

FDA CE Foptima

The PCR plate series products of Fosun Life Science are made of high-quality medical grade polypropylene, with full specifications, strong adaptability, good sealing.



FTLSLEACMBC-002

Specification	Character
0.1mL, 0.2mL	PCR tube, PCR 8-strip tube, PCR plate (non-skirt, half-skirt, full-skirt)

### Sample Tube/ Sample Container

FDA CE Foptima

Sample collection system, including varieties of containers, tubes and traps, is intended to be used in clinical lab according to the instructions provided for collection, preservation, transport and storage of dry or wet body fluids specimens, which would alleviate mess



FTLSLEACMBC-003

Specification	Character
0.5mL, 2mL, 5mL, 30mL, 40mL, 60mL, 90mL, 120mL	Bulk pack, individual pack, sterile, non-sterile

### Centrifuge Tube

FDA CE Foptima

Centrifuge tubes are made of medical-grade transparent polypropylene (PP) and have been widely used in a variety of experimental operations to meet the separation requirements of biological samples



FTLSLEACMBC-004

Specification	Character
0.5mL, 1.5mL, 2.0mL, 5.0mL, 15mL, 50mL	Conical Bottom, flat cap, bulk pack individual pack

## Pasteur Pipette / Serological Pipette

FDA CE Foptima

Pasteur pipettes are mainly used for the suction, transfer or carrying of small amounts of liquids such as cell tests, clinical tests, cloning tests, etc.



Specification	Character
Pasteur Pipette: 1.0mL, 2.0mL, 3.0mL; Serological Pipette: 1.0mL, 2.0mL, 5.0mL, 10mL, 25mL, 50mL	Bulk pack, individual pack, sterile, non-sterile

## Cell Pusher

FDA CE Foptima

Cell pusher is designed for spreading and dispersing liquids across the surface of agar culture plates. Used for performing bacterial CFU plate counts on water, milk, urine and other liquid samples.



Specification	Character
T-shaped, L-shaped, Triangle shaped	Bulk pack, sterile

## Cell Biology Consumables

### Cryogenic Tube

FDA CE Foptima

The cryogenic tube is made of transparent polymer polypropylene (PP), manufactured by a special process, so that it can withstand ultra-low temperature and ensure no leakage, and is suitable for long-term low-temperature storage of cells and tissues.



Specification	Character
1.0mL, 1.8mL, 4.5mL	Bulk pack, rack pack, sterile, non-sterile

### Cell Culture Dish/ Plate/Bottle

FDA CE Foptima

Cell culture consumables refer to the various materials and supplies used in the process of growing cells in a laboratory setting.



Specification	Character
Dish: 35mm, 60mm, 100mm, 150mm; Plate: 6/12/24/48/96 well; Flask: T25, T75, T175, T225	Bulk pack, individual pack sterile, non-sterile, TC-treat, non-TC treat

## Microbiology Consumables

### Petri Dish

FDA CE Foptima

Petri dish is made of high polymer material polystyrene. Mainly used for microbiology experiments such as laboratory inoculation, streaking, and isolation of bacteria.



Specification	Character
35mm, 60mm, 70mm, 90mm; 2 / 3 / 4 Compartments	Bulk pack, sterile

### Inoculating Loop

FDA CE Foptima

Inoculating loop is used to pick colonies during bacterial culture and perform inoculation operations. It is widely used in many disciplines such as microbial detection, cell microbes, and molecular biology.



Specification	Character
1ul, 10ul	Bulk pack, sterile

## Pathology Consumables

### Microscope Slides/ Cover Glass

FDA CE Foptima

With a great number of microscope slide specifications, Fosun Life Science are providing different choices in quality to match up with your demands. Different standards can be incorporated into production according to your requireme.



Character
Slides: 75x25mm; Cover: 18x18mm, 20x20mm, 22x22mm, 24x24mm, 22x32mm, 24x32mm, 22x40mm, 24x40mm, 22x50mm, 24x50mm, 22x60mm, 24x60mm

### Embedding Cassette

FDA CE Foptima

Fosun Life Science produce many kinds of embedding cassettes. Most models can be used with automated labeling machines



Specification	Character
Square holes, Rectangle holes, Round holes	One-piece design, Two-piece design; bulk pack

## Vacuum Blood Tube

CE

Disposable Blood Collection System offers a complete system of vacuum or non-vacuum PET/Glass/PP blood collection tubes vials containers or bottles, needles, needle holders, and safety devices also includes a complete range of products for pediatric

PlainTube; Procoagulation Tube; Gel&Clot Activator Tube; Gel&Heparin Lithium tube; Glucose Tube, Sodium Fluoride+Potassium; PTTube, 3.2%or3.8% Sodium Citrate; HeparinTube, HeparinLithiumorHeparin Sodium; EDTATube, K2EDTA, K3EDTA; ESR-Tube, 3.2%or3.8% Sodium Citrate; Oxalate;




Character
1-5mL, 6-7mL, 8-10mL

## Vaccine / Syringe & Needle / Oral & Enteral Diabetes Care / Blood Collection / Infusion Therapy

### 1ml/3ml/5ml/10ml/20ml/30ml/ 50ml/60ml LL/LS Syringe

FDA CE TGA MDSAP Fosafe

ISO 80369-7 Luer tip, compatible with standard needle  
Highly transparent barrel, graduation extended beyond




Classification
general medical masks, surgical masks medical protective masks

### Disposable Needle

FDA CE TGA MDSAP Fosafe

Made of stainless chromium nickel steel  
Bevel length according to indication



size	color
S, M, L, XL, XXL	White

### Safety Needle

FDA CE TGA MDSAP Fosafe

Designed to help prevent needlestick injuries  
Simple one-handed safety shield mechanism




size	color
S, M, L, XL, XXL	Blue, yellow, green etc

### Oral & Enteral Syringe

FDA CE TGA MDSAP Fosafe

Special tip design will not accept luer connectors, preventing medication errors and the accidental intravenous injections of oral liquids



size	color
S, M, L, XL, XXL	Blue

### 0.3ml/0.5ml/1ml Insulin Syringe

FDA CE TGA MDSAP Fosafe

Insulin syringe with attached needle for subcutaneous injection  
Highly transparent barrel with high contrast black graduation




size	color
S, M, L, XL, XXL	Blue, yellow, green etc

### Blood Collection Needle

FDA CE TGA MDSAP Fosafe

Sterile individual pack, single use only  
Comfortable double wings design, easy operation  
Safety cover reduces accidental needlesticks



size	color
S, M, L, XL, XXL	Blue

## Infusion Set

FTLSLEACVOS-007

Infusion Set  
Transfusion Set  
Extension Set  
Infusion Set With Burette



## First Aid Kit

### First Aid Box

CE

Respiratory protection: CPR face mask  
Wound dressing: Medical tape, Adhesive plaster, Dressing bandages, Burning dressing, Elastic bandages, Triangular bandage

FTLSLEACFAK-001



#### Accessory

Exam gloves, Emergency blanket, Stainless steel scissors, First aid manual

### First Aid Kit

CE

Wound dressing: Strip, Adhesive plaster, Gauze pad / Eye pad, Wound dressing, Triangular bandage, Elastic bandage, Medical tape, Cotton tip applicator

FTLSLEACFAK-002



#### Respiratory protection

CPR face mask

#### Disinfection

Alcohol pad, Povidone-iodine pad

## Medical Dressing

### Gauze Balls

FDA CE

Size: 6x6cm, 8x8cm, 10x10cm, 20x20cm, 25x26cm, 30x30cm

FTLSLEACMD-001



#### Material

17th, 20th, 24th

#### Non-sterile package

50/100/200pcs

### Non-Woven Swabs

FDA CE

Size: 5x5cm (2x2), 7.5x7.5cm (3x3), 10x10cm (4x4)

FTLSLEACMD-002



#### Weight

30g, 40g, 50g

#### Non-sterile package

100/200pcs

### Dressing Kit

FDA CE

Products list: Gauze swabs, Gauze balls, Non woven swabs, Drapes, Tweezers, Cotton balls

FTLSLEACMD-003



### Instant Cold Pack

CE

Spec: 21.5x15.3 PE+PA  
Weight: 225g  
Packing: 50

FTLSLEACMD-004



#### Ctn Size

46x30.5x23cm

#### G.W

11.3kg

## Gauze Sponges

FDA CE

Size: 5x5cm, 7.5x7.5cm, 10x10cm

FTLSLEACMD-005



#### Ply

8/12/16

#### Mesh

26/18, 17th

## Medical Tape

FDA CE

FTPPE-007



#### Size

1.25x9cm, 2.5x9cm, 5x9cm

#### Material

PE/Non woven/ZNO/Silk

## Examination/Isolation Drapes

FDA CE

Size: 20-400x20-400cm

FTLSLEACMD-009



#### Material

PE/PVC, SAP

## Surgical Operation Pack

FDA CE

Products list: Instrument table cloth, Instrument disc sleeve, Hand towel, Tap, Surgical gown, Drapes, Medical gloves, Wrap

FTLSLEACMD-006



## Gauze Bandages

FDA CE

Size: 6x600cm, 8x600cm, 10x600cm

FTLSLEACMD-008



#### Mesh

26/18, 17th

## Saline Solution

FDA CE

For adults and children of all ages. For infants under 6 months, please consult a pediatrician before use. If pregnant or breastfeeding, please consult a doctor before use.

FTLSLEACMD-010



#### size

5ml 10ml 15ml 20ml

#### Active Ingredient

SEAWATER

# Family Consumer



## Family Consumer

### Home Care

#### Multi-function Monitoring Systems

CE

measurement range: 0.1—8.0mmol/L  
Coding: Auto-coding  
Test time: 10s  
Sample type: 0.8 μL  
Hematocrit range: 20-65%  
Strip Shelf-life: 24months



Test time	Sample type
5s	0.5 μL

#### Thermometer/ Infrared Forehead & Ear Thermometer

CE

LCD size: 21x9mm  
Measurement range: 32°-42.9°  
Measurement accuracy: ±0.1°/0.2°  
Response time: 60s  
Battery: One 1.5V battery



size	color
126x18x10mm	White

#### Oximeter

CE

Product Model: 50DL  
LED display  
Display of SpO<sub>2</sub>, PR and bar graph  
Low battery indication

Product Model: 50DL  
Auto power on color display  
Display of SpO<sub>2</sub>, PR, bar graph and pulse waveform



Power supply	Continous working time	Packing	CTN Size
two batteries	30h	200 pcs/box	433x425x353mm 10kg

#### Arm Blood Pressure Monitor

CE

Product Model: C02  
Display: LED screen  
ABS shell  
Carry case of Cuff  
Power: Internal-DC 4.5V (3\*AA)

Product Model: 65B  
Display: Large LCD screen  
Voice broadcast  
Cuff Size: 23-33cm  
Power: DC 6V/500mA



Display	Memory	Momery	Packing
LED screen	2x120 Memories	2x192 sets	30 pcs/ctn, 20 pcs/ctn

#### Blood Glucose Monitor

CE

Enzyme: glucose oxidase  
Sample size: about 0.5 μL  
Testing time: about 5 seconds  
Hematocrit range: 20-70%  
Battery: 250mAh 3.7V rechargeable battery

Enzyme: glucose dehydrogenase  
Sample size: about 0.6 μL  
Testing time: about 5 seconds  
Hematocrit range: 30-55%  
Operation temperature: 5°C - 45°C (41°F - 113°F)



Operation temperature	Memory
5°C - 45°C (41°F - 113°F)	500 results with time and date

# Family Consumer

## ■ Oxygen Concentrator

### 3L/5L/10L Oxygen Concentrator

FDA CE

Medical standard, high oxygen concentration  
365 days \* 24 hours constant oxygen supply  
Oxygen concentration and nebulization  
Real-time oxygen concentration & flow monitoring system  
Microcomputer (MCU) intelligent control



## ■ UV-Ozone Disinfection Cart

### UV Disinfection Cart

CE

Medical UV Disinfection Cart  
Application Area: ≤20 m<sup>2</sup>  
Disinfection Timing: 0~120mins  
Features: It utilizes ultraviolet to eliminate bacteria fungi, mites, and coronavirus



FTFCUVDC-001

### UV-Ozone Disinfection Car

CE

Medical UV Disinfection Cart  
Application Area: ≤50 m<sup>2</sup>  
Disinfection Timing: 0~150mins  
Features: Ultraviolet plus ozone greatly improves the disinfection efficiency



FTFCUVDC-002

Model	Packing
SX	1040*210*270mm

Model	Packing
SX-150	1040*210*270mm

# FOSUN TRADE

Provide Global Customers with  
One-stop Integrated Trade Solution



Bund Finance Center,600 Zhongshan No.2 Road(E),Shanghai 200010,P.R.China

Tel +86 21 23156666

Fax +86 21 61303288

www.fosuntrade.com