

Nasogastric tube is a single lumen, small bore feeding tube. It is more appropriate for administration of medication or nutrition.



Product Code	Size	Length (mm)	Color Code	Certificate
NA04P	4Fr	400	Red	CE, ISO, CFDA, FSC
NA05P	5Fr		Grey	
NA06P	6Fr		Light Green	
NA08P	8Fr		Light Blue	
NA10P	10Fr	1200	Black	
NA12P	12Fr		White	
NA14P	14Fr		Green	
NA16P	16Fr		Orange	
NA18P	18Fr		Red	
NA20P	20Fr		Yellow	
NA22P	22Fr		Purple	



The ENFit connector features a unique design to reduce the risk of enteral feeding tube misconnections and improve patient safety. It addresses "patient side" connections between feeding tubes, administration sets, medication, flush and bolus feeding syringes and other enteral devices, but does not allow connectivity with any other connector for any other clinical use.





The anterior section of the guide wire is passivated, which will reduce the risk of causing damage to the esophagus and gastric mucosa. Stainless steel material is convenient for medical staff to use; multi-purpose connector can be connected with different tubes.

