


UNIVERSAL SUTURES
"HEALING BEYOND COMFORT"




Sutures
Mesh
Bone Wax

AN ISO & CE CERTIFIED SURGICAL SUTURES AND MESH MANUFACTURING COMPANY
We Manufacture complete range of Surgical Sutures, Hernia Mesh & Bonewax

Needle Shape



Needle Type:




International Coding Pattern

Code 1	Product	Code 2	USP Size	Code 3	Code 4	Needle Curvature	Code 5	Needle Type	Code 6
1	Braided & Coated Polyglycolic Acid Undyed	A	3	SUTURE LENGTH in cm	C	3/8 Circle	Z	Round Body / Taper Point	NEEDLE LENGTH in mm
2	Braided & Coated Polyglycolic Acid	B	2		D	1/2 Circle	Y	Cutting	
3	Monofilament Polyamide (Nylon)	C	1		A	Straight	X	Reverse Cutting	
4	Catgut Chromic	D	0		E	5/8 Circle	K	Taper Cutting	
5	Black Braided Silk	E	2-0					T	Trocar Point
6	Braided Polyester	F	3-0					R	Blunt Point
7	Catgut Plain	G	4-0					S	Spatula
8	Monofilament Polypropylene	H	5-0						
9	Monofilament Polydioxanone	I	6-0						
10	Braided & Coated Polyglactin 910	J	7-0						
11	Monofilament Polyglecaprone 25	K	8-0						
12	Braided & Coated Polyglactin 910 Undyed								

Example:

Braided & Coated Polyglycolic Acid	Size 1	75cm	1/2 Circle	Round Body	40mm
------------------------------------	--------	------	------------	------------	------



UNIGLYDE FAST

Braided & Coated Polyglycolic Acid (Fast Absorbable)
Absorbable Surgical Suture U.S.P.



UNIGLYDE FAST - Braided & Coated Polyglycolic Acid Suture is an absorbable, sterile, synthetic surgical suture, composed of low molecular weight homopolymers of Glycolide (100%).

Braided & Coated Polyglycolic Acid Suture meets all the requirements as per United States Pharmacopeia for Absorbable Surgical suture (Synthetic)

In-house needle manufacturing provides you the better needles with optimum sharpness which finally delivers smooth riding on tissues during suturing.

Indications of Polyglycolic Acid Suture

Polyglycolic Acid (PGA Suture) suture is indicated for use in general soft tissue approximation and/or ligation, including use in ophthalmic procedures, but not for use in cardiovascular or neurological tissues.

- Smooth Surface provides no tissue trauma
- Secure Knot holding capacity
- Best suture for critical wound support
- Provides Best tissue support even after smaller diameter
- Hydrolytic absorption ensures no tissue reaction

Additional Information:

Characteristics	:	Synthetic Absorbable Surgical Suture
Type	:	Braided
Material	:	100% Polyglycolic Acid synthetic absorbable braided surgical suture
Coating	:	Polycaprolactone & calcium stearate
Color	:	Undyed
Support	:	7 to 10 days
Absorption	:	Absorption takes place through hydrolysis completely in 40 to 42 days.
Strength	:	75% tensile strength is retained for 3 days and 50% till 5 days.
USP Range	:	6/0 – 2
Sterilization	:	EO (Ethylene Oxide)
Shelf Life	:	3 Years
Product Characteristics	:	Excellent tensile strength, Excellent knot tying, Minimum tissue reaction
Packaging	:	In aluminium foil packages with or without needle
Needle Type	:	Round bodied, cutting edge, taper point, straight, blunt point
MFG. LIC No	:	MFG/MD/2020/000109
OEM. LIC No	:	N. Code: KA/DEVICE/MFG/MD/2020/000109

UNIGLYDE FAST

Braided & Coated Polyglycolic Acid (**Fast Absorbable**)
Absorbable Surgical Suture U.S.P.



Code No.	SUTURE		NEEDLE DESCRIPTION
	USP SIZE	LENGTH	
1H75DZ16 Lot 47	5-0	75 cm	1/2 Circle Round Body 16cm
1H75DZ20 Lot 48	5-0	75 cm	1/2 Circle Round Body 20cm
1G75DX16	4-0	75 cm	1/2 Circle Reverse Cutting 16cm
1G75DZ16 Lot 43	4-0	75 cm	1/2 Circle Round Body 16cm
1G75DZ19	4-0	75 cm	3/8 Circle Reverse Cutting 19cm
1G75DZ20 Lot 44	4-0	75 cm	1/2 Circle Round Body 20cm
1F75DZ16 Lot 15	3-0	75 cm	1/2 Circle Round Body 16cm
1F120DZ17 Lot 196	3-0	120 cm	1/2 Circle Round Body 17cm
1F75DZ17 Lot 207	3-0	75 cm	1/2 Circle Round Body 17cm
1F70DZ17	3-0	70 cm	1/2 Circle Round Body 17cm
1F90DZ20	3-0	90 cm	1/2 Circle Round Body 20cm
1F75DZ20 Lot 17	3-0	75 cm	1/2 Circle Round Body 20cm
1F75DZ26 Lot 198	3-0	75 cm	1/2 Circle Round Body 26cm
1E75DZ20 Lot 34	2-0	75 cm	1/2 Circle Round Body 20cm
1E75DZ26 Lot 35, 197	2-0	75 cm	1/2 Circle Round Body 26cm
1E75DZ30 Lot 205	2-0	75 cm	1/2 Circle Round Body 30cm
1E75DZ36 Lot 200	2-0	75 cm	1/2 Circle Round Body 36cm
1D75DZ26 Lot 28	0	75 cm	1/2 Circle Round Body 26cm
1D70DZ26	0	70 cm	1/2 Circle Round Body 26cm
1D90DZ30 Lot 201	0	90 cm	1/2 Circle Round Body 30cm
1D90DZ36 Lot 202	0	90 cm	1/2 Circle Round Body 36cm
1D90DZ40 Lot 30	0	90 cm	1/2 Circle Round Body 40cm
1C90DZ30 Lot 204	1	90 cm	1/2 Circle Round Body 30cm
1C90DZ36 Lot 49	1	90 cm	1/2 Circle Round Body 36cm
1C75DZ36 Lot 199	1	75 cm	1/2 Circle Round Body 36cm
1C90DZ48 Lot 203	1	90 cm	1/2 Circle Round Body 48cm
1B100DZ48 Lot 206	2	100 cm	1/2 Circle Round Body 48cm