


Redon Drainage Tube

- Provided with **radio-opaque line** and smooth eyes
- Facilitates **efficient drainage with minimal tissue trauma**
- Cross perforated length of tube: 14 cm
- **Markings on the tube at 19 cm from the distal tip for monitoring insertion depth and drainage**
- Soft, frosted & kink resistant **PVC** tubing
- Length: 50, **70**, 120 & 170 cm
- Sizes: 8, 10, 12, **14**, **16** & **18** FG
- **DEHP Free**



Redon tube dimensional parameters are as per below table:

Size	O.D. ± 0.10 (mm)	I.D. ± 0.10 (mm)	Perforated Length (cm)		Number of perforations at equal distance	
			50cm Tube	70cm Tube	50 cm Tube	70cm Tube
FG 8	2.70 mm	1.60 mm	14 cm	14 cm	24 pairs	24 pairs
FG 10	3.33 mm	2.10 mm	14 cm	14 cm	24 pairs	24 pairs
FG 12	4.00 mm	2.50 mm	14 cm	14 cm	24 pairs	24 pairs
FG 14	4.70 mm	2.90 mm	14 cm	14 cm	24 pairs	24 pairs
FG 16	5.33 mm	3.60 mm	14 cm	14 cm	26 pairs	26 pairs
FG 18	6.00 mm	4.20 mm	14 cm	14 cm	26 pairs	26 pairs

Sterilization Method:

Sterilized using Ethylene Oxide.

Packaging Characteristics:

Individual unit pouch.