


TAIZHOU KANGJIAN MEDICAL EQUIPMENT CO.,LTD

TAIZHOU KANGJIAN MEDICAL EQUIPMENT CO.,LTD.TIANCHANG BRANCH



台州康健医用器械有限公司
Taizhou kangjian Medical Equipment Co.,Ltd.

地址：浙江省玉环市机电工业园区
Address: The Machine Electricity Zone of Yuhuan City,
Zhejiang Province 317600, P.R.China
TEL : 0086-576-87299799 87225930
FAX : 0086-576-87239799 87239885
P.C : 317600
Http://www.kangjiancn.com
E-mail: webmaster@kangjiancn.com

台州康健医用器械天长分公司
Taizhou kangjian Medical Equipment Co.,Ltd.

Tianchang Branch
地址：安徽省天长市铜城镇
Address: Tongcheng Town, Tianchang City, Anhui, China
TEL : 0086-550-2397575
FAX : 0086-550-2397677

KANGJIAN

Taizhou kangjian Medical Equipment Co.,Ltd.
台州康健医用器械有限公司

台州康健医用器械有限公司
Taizhou kangjian Medical Equipment Co.,Ltd.
台州康健医用器械天长分公司
Taizhou kangjian Medical Equipment Co.,Ltd.Tianchang Branch



企 業 簡 介

台州康健医疗器械有限公司，现坐落于东海之滨—浙江省玉环市机电工业园区，创建于1993年。企业占地面积20000M²，建筑面积15000 M²，库房面积8000 M²，其中符合GMP十万级净化车间面积8000 M²。现有固定资产2500万元。职工300人，其中工程技术人员30名，质量管理体系内审员5名。

公司主要从事一、二类医疗器械的生产与销售。主要产品有：一次性使用无菌阴道扩张器、吸烟式扩阴器、带光源扩阴器、脚踏式扩阴器、一次性阴道冲洗器、带光源肛镜、面罩系列、鼻氧管、一次性使用妇科包、一次性使用宫颈取样器、透气胶贴、棉签、碘伏棉签、检查手套、一次性使用薄膜手套、一次性使用压舌板、扩阴器等。

公司生产、检测设备和仪器先进、齐全，完全满足生产及检测的需要。从事产品出厂检验人员全部经过国家药品监督管理局所指定的监督检验中心培训，考试合格，持证上岗。

公司秉承“依法管理、质量第一、信誉至上”的质量方针，以人为本，建立健全并保持质量管理体系的有效性。质量管理体系及产品已获得ISO13485、ISO9001、CE及FDA证书，确保产品安全、有效。过硬的产品质量、良好的企业信誉，赢得了丰厚的市场回报。产品已畅销全国，并出口欧洲、中东、东南亚、澳洲、非洲、美洲等国家和地区。

热忱欢迎业内各界朋友惠顾指导，共谋发展，为人类的健康事业贡献我们的力量。



Company Brief Introduction

Taizhou Kangjian Medical Equipment Co., Ltd. was established in 1993, located in Yuhuan, Zhejiang province, China.


We are the professional enterprise for manufacturing various disposable medical devices. The main products including Disposable Vaginal Speculum, Vaginal Speculum with Smoke Extractor, Vaginal Speculum with Light Source, Self-destruct Vaginal Speculum, Gynecological sets, Cervical Sampler, Vaginal Irrigator, Oxygen Mask Series, Nasal Oxygen Cannula, Anoscope with Light Source, Nostril dilator, Tongue Depressor, Medical Cotton-tipped Applicator Stick, Iodine Cotton Swab, Iodine Cotton Ball, Sterile fixed sticking cloth for transfusion needle, PE glove, etc.

The products have been well sold to all of the country, and also exported to Europe, Asia, America, South America middle east, Australia, Africa, etc.

We have established and maintained the quality management system, and have achieved the certificates of the ISO13485, CE and FDA.

Welcome friends from all circles to patronize us for common development and contribute our honesty and affection to the healthy course of human being.






Cervical Brush (Item no.: KJ-CB02)

Features:

- * Length:202mm / Brush width:17.5mm / Handle width:5.1mm / Brush height:25.5mm
 - * Material: polyethylene
 - * Blue handle or any other colors
 - * Intended use: used in collecting cervical exfoliated cells for smear test.
 - * Sterile EO
 - * Packaging: 1 piece/bag,100pcs/box,2000pcs/outer carton
- Inner box size: 265x130x90mm
Outer carton size: 680x300x390mm



Cyto Brush (Item no.: KJ-CB03)

Features:

- * Length:170mm / Brush width:4.1mm / Handle width:5.0mm / Brush height:20.7mm
 - * Material: polyethylene handle & Nylon bristles with 304 stainless steel
 - * Blue handle or any other colors
 - * Intended use: used in collecting cervical exfoliated cells for smear test.
 - * Sterile EO
 - * Packaging: 1 piece/bag,100pcs/box,2000pcs/outer carton
- Inner box size: 265x120x80mm
Outer carton size: 620x300x350mm



Gynecological Kit (Item no.: KJ-GK02)

Features:

Including items:

- I. Cusco type vaginal speculum side screw type - 1 piece
- II. disposable latex gloves - 2 pcs.
- III. Cervex-brush combi - 1 pc.
- IV. sheet (bedding) with dimensions 40cmx60cm - 1 pc
- V 1.individually packed (set)

* Intended use: used in checking gynecology disease in hospitals or clinics for smear test. Before use, the user must read the instructions for use.

* Sterile EO

* Packaging: 1 piece/bag,100pcs/carton


Outer carton size: 560*380*350mm



Disposable Vaginal Speculum (Item No.: KJ-V02)

Features:

- * Material: Medical grade transparent Polystyrene
- * Structure: upper duckbill & down duckbill & nut & screw
- * Intended use: used in checking and treatment gynecology disease in hospitals or Clinics.
- * Sterile, Non-toxic, Pyrogen free, Hypoallergenic
- * Packaging: 1 piece/bag, 100pcs/carton
Outer carton size: 500*380*280mm



3.1 Vaginal Speculum Lateral Screw Type



L(white nut) M(red nut) S(yellow nut)

Basic Size:

Specification	Basic Size (mm)			
	l	w	b	h
Large	120±2	34±1	32±1	57±3
Medium	110±2	34±1	28±1	55±3
Small	104±2	34±1	20±1	53±3