



Histo-Line Laboratories (/en)

Code : **LP0001**



ASSISTANCE?

Request information ([mailto:customer-care@histoline.com?Subject=Request information - product:Super Pap Pen](mailto:customer-care@histoline.com?Subject=Request%20information%20-%20product:Super%20Pap%20Pen))

T +39 02 55230061 (tel:+390255230061)

| Description | Packaging | Datasheet & SDS |
|---|-----------|-----------------|
| <p>The Super PAP Pen is useful for immunohistochemical and fluorescent staining methods.</p> <p>It prevents the waste of valuable reagents by creating a water repellent circle around the section.</p> <p>In addition, it is now chemically formulated to withstand the rehydration steps performed in alcohol, as well as the high temperatures required for denaturation during in situ hybridization.</p> | | |