



Double-pedal footswitch for ARC, ARC 250 / 303, cable 4 m

BOWA ARC foot switch series for ARC, ARC 303, ARC 250

Scope of delivery

Incl. instructions for use

Product number: 901-022

Unit: 1 piece

Type: with lateral stop

 [Reset selection](#)

Inquiry

Description Specification Download

ARC footswitch without switch

The BOWA foot switch series is available in different models depending on the ARC device system used. Depending on the foot switch type chosen, a cut and coag pedal and a cable length of 4 m are available for a generous action radius. Effective cleaning and safety are ensured with the ability to immerse the device completely and the explosion protection.

Scope of delivery

Incl. instructions for use

Do you still have any questions regarding the product?
We will be happy to help you!

Get inquiry now

Have you seen?



Electrosurgical unit ARC 303

Powerful top-class electrosurgical cutting and ...

1 piece



Electrosurgical unit ARC 250

Top quality cutting and coagulation in electros...

1 piece

