HIGH VACUUM, HIGH FLOW HOSPITAL SUCTION ASPIRATORS









Customized Hospi Plus MPR with 1x2 I jar, 1x5 I jar and safety trap Multipurpose stainless steel rail for maximum versatility allows to

hang any accessory. Suitable for MPR models only 28154. 28155. 28156

Rings for: 2 or 5 l autoclavable jars,

2 or 3 I Flovac

Designed for professional aspiration of bodily fluids, tissues or bones of patients during or after surgery.

Oiless and maintenance free piston type pumps provide high performances. Excellent suction capacities and max vacuum built up within a few seconds. Available with vacuum gauge and autoclavable jars made of Makrolon® with 200 ml (2 l jar) or 500 ml (4 or 5 l jars) graduation optional 2-3 litres Flovac jars with disposable liners, not compatible with safety trap bottle. Wide range of versions with different features

- 60 l/m (Clinic Plus) or 90 l/m (Hospi Plus) flow rate;

- 2, 4 or 5 l jars;

Power Cord with Schuko plug

Multilingual manual: GB, FR, IT, ES, DE, GR, PL, RO

and HOSPI PLUS

Membrane keyboard for 28119-20, 28160,

28162, 28195

Clinic Plus Full and Hospi Plus Full with foot switch and electronic flow direction regulator TECHNICAL SPECIFICATIONS

continuous 90 l/min

- MPR system for versatility;

- footswitch and flow direction diverter.

High standards of safety in overflow protection system Double security valve integrated in the jar and hydrophobic filter (all models), safety trap bottle (only MPR models). MPR (Multi Purpose Rail) models

This system enhances the versatility for easy and quick exchange of different accessories, with no need for tools. Equipped with five connections for rings of various diameters to fit jars of different sizes and types (2I, 3I, 5I), cannula holders and safety trap bottle.



Multipurpose rail with 5 locking slots offers versatility to operators in changing accessories

Safety trap bottle

	Voltage		SUCTION ASPIRATORS	Autocia-	Foot	FIOW
	230 V	110 V	CLINIC and HOSPI PLUS	vable Jars	switch	direction diverter
	28194	28118	CLINIC Basic	2x2 l	-	-
ì	28195	-	CLINIC Full	2x2 l	yes	yes
ı	28196	-	CLINIC Basic	2x4 l	-	-
ı	28198	28119**	CLINIC Full	2x4 l	yes	yes
ı	28159	-	CLINIC Basic MPR	2x5 l	-	-
ı	28160	-	CLINIC Full MPR	2x5 l	yes	yes
ı	28200	-	HOSPI Basic	2x2 l	-	-
ı	28203	-	HOSPI Full	2x2 l	yes	yes
ı	28201	-	HOSPI Basic	2x4 l	-	-
ı	28204	28120**	HOSPI Full	2x4 l	yes	yes
ı	28161	-	HOSPI Basic MPR	2x5 l	-	-
ı	28162	-	HOSPI Full MPR	2x5 l	yes	yes

Operating voltage:	← 230 V - 50/60 Hz or 110 V - 6
Maximum suction:	- 0.90 bar (675 mmHg)
Operating cycle:	continuous
Flow:	60 l/min
Power:	230 VA

CLINIC PLUS

Clinic Plus Basic and Hospi Plus Basic without

foot switch and electronic flow direction regulator

385 VA Noise level: 51.7 dBA 46.4 dBA Norms: EN 60601-1-1, EN 60601-1-2 93/42 EEC

HOSPI PLUS MAIN APPLICATIONS EMERGENCY DEPT. **GENERAL SURGERY** GYNAECOLOGY NEUROSURGERY **OBSTETRICS OPERATING THEATRE**

DENTAL PRACTICE

STANDARD ACCESSORIES CLIN		Clinic,	Clinic	Clinic	Hospi basic MPR	Hospi full MPR
121°C autoclavable jar with overflow valve system (2, 4 or 5 l)	2	2	2	2	2	2
Rings to accommodate 5 l jars	-	-	2	2	2	2
Safety Trap Bottle (220 ml)	-	-	1	2	1	2
Antibacterial & Hydrophobic Filter (single-patient)	1	1	1	1	2	2
Autoclavable silicone tubes Ø 8x14 mm length 150 cm	1	1	1	1	2	2
Conical Connector Ø 10-11-12mm	1	1	1	1	2	2
Air suction inlet	1	1	1	1	2	2
Footswitch with intermittent or continuous operation	-	1	-	1	-	1
Change-Over System from jar to jar by soft-touch keys	-	1	-	1	-	1

20204	20120	HOSFITUII	2A4 I	y C3	yes			
28161	-	HOSPI Basic MPR	2x5 l	-	-			
28162	-	HOSPI Full MPR	2x5 l	yes	yes			
Allows to direct suctioned liquids to any of the 2 jars **with 2x5 l jars								
GIMA code	ACCESSORIES FOR CLINIC/HOSPI PLUS MPR							
28150	Safety trap bottle 220 ml with tube - spare							
28151	Multipurpose stainless steel rail - optional							
	Ring for 2 I autoclavable jar - optional							
28155	Ring 1	for 5 I autoclavable	jar - spa	are				

28156 Ring for 2 or 3 | Flovac jar - optional **MEMBRANE KEYBOARD - spare**

28170 Mem. Keyboard for 28198 **28171** Mem. Keyboard for 28119-20, 28160, 28162, 28195 28172 Mem. Keyboard for 28203, 28204

JAR HOLDERS - spare **28108** Jar holder 2 l* for 28194-5, 28118, 28200, 28203

28109 Ring adapter for 28108 *Needs 28109