

### (1) Multi Functions

- ClearCoajet is a combination of Injection catheter and Hemostat.
- Marking, Injection, Coagulation are all possible with this device for EMR, ESD procedure.

### (2) Back-cut Needle design

• It has advantages to Insertion&injection saline because the shape of Needle is sharp and I.D(Inner Diameter) is large.




Injection without



### SPECIFICATION

Model	Needle O.D	Needle Thickness	Needle Length	Sheath 0.D	Working Length	Number of Use	Validity of Sterilzation
FM-CJ0002(1800) FM-CJ0002(2300)	23G	0.1	4mm	¢2.4mm	1,800mm 2,300mm	Single Use Only	3 years

## ClearHemograsper



### (1) Rotation with the handle

- Enhanced the Rotation, and the users can perform Hemostasis for bleeding at demanding location more easily.
- Jaw and Sheath are rotated separately.

### (2) Installed irrigation port

• Provide convenience to users that use the scope without Water Jet function.

# Irrigation Function Syringe or water pump irrigation port

### SPECIFICATION

Model	Opening width	Sheath O.D	Type of Rotation	Working Length	Number of Use	Validity of Sterilzation
FM-HG0001(1650) FM-HG0001(1820) FM-HG0001(2000)	5mm	¢2.3mm	Rotational	1,650mm 1,820mm 2,000mm	Single Use Only	3 years

Imagine What a Woderful difference You have 7