

DAHLHAUSEN®



Bacterial/Viral Filter, with HME

Bacterial/Viral Filter with HME and CO₂-port. The purpose of the filter is to protect the patient and to humidify the inhaled air. The filter is placed between the endotracheal tube/tracheal cannula and the Y-piece of the breathing system.

Shape (Ventilation Filter): Round

ISO style (patient end): Male/Female

Inner diameter (patient end): 20/15 mm

Sterile: No

Art.no.	Application area	Particle retention	ISO style (machine end)	Inner diameter (machine end)	Max. tidal volume (Breathing Filter)
4300906120	Child	99,99 %	Male	13 mm	900 ml
4300906020	Adult	99,999 %	Male/Female	22/13 mm	1.500 ml
4300905708	Adult	99,999 %	Male/Female	22/13 mm	1.500 ml

P. J. Dahlhausen & Co. GmbH
Emil-Hoffmann-Str. 53
50996 Köln