

| Ref. no. | Product | Units | Illustration | | | |
|--|--------------------------|-------|--------------|--|--|--|
| - Fibrin Set (Preparation Kit and Application Kit) | | | | | | |
| VS 302 | Fibrin Set | 10 | + | | | |
| VS 312 | Fibrin Set - Concorde | 10 | + | | | |
| VS 322 | Fibrin Set - Co-Delivery | 10 | + | | | |
| VS 323 | Fibrin Set - Endoscopic | 10 | + | | | |
| - PRF® Set (Preparation Kit and Application Kit) | | | | | | |
| VS 400 | PRF® Set | 10 | + | | | |
| VS 410 | PRF® Set - Concorde | 10 | + | | | |
| VS 420 | PRF® Set - Endoscopic | 10 | + | | | |
| VS 422 | PRF® Set - Co-Delivery | 10 | + | | | |
| - Other | | | | | | |
| VS 510 | Vivostat® Split Kit | 20 | | | | |



| Ref. no. | Product | Units | Illustration | | | | |
|--|--|-------|--------------|--|--|--|--|
| - Application Kit (Including all necessary components for Application) | | | | | | | |
| VS 305 | Spraypen® Kit | 4 | 200,000 | | | | |
| VS 315 | Spraypen® Kit - Concorde | 4 | | | | | |
| VS 325 | Endoscopic Kit ¹ | 4 | | | | | |
| VS 335 | Spraypen® Kit - Co-Delivery | 4 | Name . | | | | |
| VS 345 | Endoscopic Kit - Straight ² | 4 | | | | | |
| VS 355 | Endoscopic Kit - Co-Delivery ¹ | 4 | | | | | |
| - Preparat | - Preparation Kit (Including all necessary components for Preparation) | | | | | | |
| VS 306 | Fibrin Preparation Kit | 10 | | | | | |
| VS 406 | PRF® Preparation Kit | 10 | | | | | |
| - Durables | | | | | | | |
| APL 400 | Applicator Unit | 1 | 1 | | | | |
| APL 404 | Applicator Unit Co-Delivery ³ | 1 | T ce ec | | | | |
| PRO 800 | Processor Unit | 1 | (v | | | | |

 $^{^1}$ Only to be used together with VS 220 - Endoscopic Applicator Handle 2 Can be used together with VS 220 and scope solutions such as Bronchoscope, Arthroscope, Laparoscope etc. 3 All applications need to be controlled with the Vivostat $^{\tiny @}$ Foot Switch, VS 222

Vivostat® Product List 2020



| Ref. no. | Product | Units | Illustration | | | |
|--------------------------------------|------------------------------|-------|--------------|--|--|--|
| VS 220 | Endoscopic Applicator Handle | 1 | | | | |
| VS 222 | Foot switch | 1 | | | | |
| - Installation and Demonstration Kit | | | | | | |
| VS 280 | Installation Kit | 1 | + | | | |
| VS 290 | Demonstration Kit | 1 | + | | | |