



SCALE/SCALA 1:1

1 URINE VACUUM TUBES PROVETTE SOTTOVUOTO PER URINE

Sterile vacuum test tubes for chemical analysis of urine (cylindrical bottom) and microscopic examination of urine sediment (conical bottom). For urine collection tube without additive, it's recommended that the analysis be carried out as soon as possible after collection.

Provette sterili sottovuoto per l'analisi chimica delle urine (fondo cilindrico) e l'esame microscopico del sedimento urinario (fondo conico). Per il prelievo di urina con provette senza additivo, è consigliabile che l'analisi sia effettuata il più presto possibile dopo la raccolta.



COD.	DIM. MM	BOTTOM SHAPE FONDO	FILLING VOLUME VOLUME DI RIEMPIMENTO	CAP COLOUR COLORE TAPPO
35150	Ø 16 x 100	Round / Cilindrico	9 ml	Yellow / Giallo
35100	Ø 16 x 100	Conical / Conico	9.5 ml	Yellow / Giallo

NEW items

2 URINE VACUUM TUBES WITH PRESERVATIVE PROVETTE SOTTOVUOTO PER URINE CON PRESERVANTE

Sterile vacuum test tubes with urine preservative for samples of urine that can not be analyzed within two hours of collection and urine samples with parameters must be stabilized that would otherwise be altered. The preservative maintains urine for up to 24/48 hours at room temperature. The preservative helps to prevent overgrowth without causing toxicity to existing pathogens.

Provette sterili sottovuoto con preservante per i campioni di urina che non possono essere analizzati entro le due ore dalla raccolta e per campioni di urina con parametri che devono essere stabilizzati al fine di non essere alterati. Il preservante mantiene le urine per un massimo di 24/48 ore a temperatura ambiente. Il preservante aiuta a prevenire la crescita eccessiva senza causare tossicità di agenti patogeni esistenti.



COD.	DIM. MM	BOTTOM SHAPE FONDO	FILLING VOLUME VOLUME DI RIEMPIMENTO	CAP COLOUR COLORE TAPPO
35250	Ø 16 x 100	Round / Cilindrico	9 ml	Yellow / Giallo
35200	Ø 16 x 100	Conical / Conico	9.5 ml	Yellow / Giallo

NEW items