

# iM80 Patient Monitor



**iM80**  
Patient Monitor



## About Edan

Edan is a healthcare company dedicated to improving the human condition around the world by delivering value-driven, innovative and high-quality medical products and services. For over 20 years, Edan has been pioneering a comprehensive line of medical solutions that address a broad range of healthcare practices including:

- Diagnostic ECG
- Patient Monitoring
- OB/GYN
- Ultrasound Imaging
- Point-of-Care Testing
- In-Vitro Diagnostics
- Veterinary

Healthcare professionals around the world depend on Edan's breakthrough medical technologies and outstanding customer support.



**U.S. and Canada inquiries:**  
EDAN Diagnostics, Inc. | 9833 Pacific Heights Blvd., Suite E/F  
San Diego, CA 92121-4707 | +1.858.750.3066  
www.edandiagnostics.com | edan-info@edandiagnostics.com

**Inquiries outside of the U.S. and Canada:**  
Edan Instruments, Inc. | No.15 Jinhui Rd., Jinsha Community,  
Kengzi Subdistrict, Pingshan District, Shenzhen | 518067 P.R. China  
+86.755.26898326 | www.edan.com.cn | info@edan.com.cn

© 2018 Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission. Not all products or features are available in all countries, contact Edan for local availability.



ENG-PM-iM80-V2.7-20180515



# iM80

## Patient Monitor



15"


Maximum 13 Waveforms



Shortcut Menu Night Mode



Pitch Tone Thermal Recorder




### Standard Parameters:

- 3/5-lead ECG, HR, RESP, EDAN SpO<sub>2</sub>, EDAN NIBP, PR, 2-TEMP

### Optional Parameters:

- 12-lead ECG, Suntech NIBP, 4-IBP(with ICP support), C.O., Respiration CO<sub>2</sub> (Mainstream and Sidestream), MASIMO AG/O<sub>2</sub>, Nellcor OxiMax™ SpO<sub>2</sub>, EDAN G2 CO<sub>2</sub> (Sidestream)

### Multiple Display Modes



### EtCO<sub>2</sub> (Mainstream/Sidestream) for Intubated/Non-Intubated Patients



#### EDAN G2 CO<sub>2</sub> (Sidestream)

- Superior water trap design for accurate monitoring
- iCARB™ algorithm with intelligent CO<sub>2</sub> pseudo wave identification technology
- Multiple sampling accessories as options for adult, pediatric and neonate patients

#### Respiration CO<sub>2</sub> (Mainstream/Sidestream)

- Plug & Play
- Dehumidification tube instead of water trap design
- No need to calibrate on a regular basis
- Low sampling rate of 50 ml/min

### Anesthetic Gas/O<sub>2</sub> (Mainstream/Sidestream)



#### MASIMO Anesthetic Gas/O<sub>2</sub> (Mainstream/Sidestream)


- Unique Plug & Play module for mainstream anesthetic gas monitoring
- The no moisture sampling line design suitable for all clinical applications
- Low sampling rate of 50ml/min for all types of patients

#### Dräger AG (Sidestream)

- Unique water cup design
- Widely Compatible sampling tubes with luer connector

### Reliable Algorithms

- iSEAP™ ECG algorithm optimized for arrhythmia detection, pacemaker detection, and HR measurement
- SEMIP® 12-lead diagnosis algorithm with 208 kinds of analysis results
- iMAT™ SpO<sub>2</sub> algorithm with outstanding motion resistance and low perfusion resistance performance
- iCUFFS™ NIBP algorithm optimized for cardiac patients, hypertensive patients, and neonatal patients



Wall Mount Rolling Stand