

# WOUND DRESSING

## FY1101 Gauze Swabs

The pads are folded all by machine, pure 100% cotton yarn ensure the product soft and adherent. Superior absorbency makes the pads perfect for absorbing blood and exudates. In accordance with the customers' requirements we can produce different kinds of pads, such as folded, and folded, with x-ray and non x-ray. The adherent pads are perfect for operation. General specification: 2"×2", 3"×3", 4"×4" PLY: 8ply, 12ply, 16ply Mesh: 19×10, 19×15, 20×12, 26×18 EO Sterile / Non-Sterile



## FY1105 Non-woven Swabs

Non-woven swabs from 70% polyethylene and 30% terylene. The pad is designed to wick away fluids and disperse them evenly. Product have been cut like "O" and "Y", so it is easy to use. It is mainly used to absorbing blood and exudates during operation and cleaning the wounds. General specification: 30g/m<sup>2</sup>, 40g/m<sup>2</sup>; 2"×2", 3"×3", 4"×4" PLY: 4ply, 8ply EO Sterile / Non-Sterile



## FY1103 Laparotomy Sponge

Laparotomy Sponge for absorbing blood and exudates or for cleaning wound. Cotton Yarn: 40's x 40's It can prevent the wound from being effect. Our laparotomy sponges are sewn and cross-stitched. General specification: 45×45cm/30×30cm, PLY: 8ply, 16ply Mesh: 19×15/26×18 EO Sterile / Non-Sterile Prewashed/ Non-Prewashed



## FY1102 Gauze Roll

Material : 100% cotton Cotton Yarn: 40's×40's Weaving according to customer requirements, high pressure, high temperature degreasing, bleaching, refined but become, with x-ray or without x-ray Products without fluorescence, soft, has the very strong absorption ability, is the first choice for medical units use products, quality in line with the British pharmacopoeia standards, are widely used in other markets General specification: 36"×100yard-4ply, 36"×100m-4ply, 36"×1000m-4ply; Mesh: 19×11/19×15

